



PRODUCT CATALOG



UCT ENVIRO

Founded in 1986, UCT has grown to be a respected leader in the drug testing, pharmaceutical, clinical, environmental and agricultural industries. Our wide range of highly reproducible solid-phase extraction columns allow the chromatographer a consistent extraction technique, and our expertise in silane manufacturing allows greater control of the chemical processes involved in producing our high quality bonded phases.

We manufacture our complete product line of bonded silica sorbents, packaged in a variety of formats, including SPE columns, 96 & 48 well plates, Universal cartridges and micro-centrifuge tubes. We also offer a variety of SPE accessories including derivatizing reagents, GC liners, and manifolds. Our commitment to ensuring the satisfaction of our customers is accomplished by delivering on our promises: top-quality, dependable solid-phase extraction and chromatography products, and unmatched technical support.

A GREENER EARTH

ORGANIZATIONS WE SUPPORT: ARBOR DAY FOUNDATION, AUDUBON SOCIETY, SIERRA CLUB

Here at UCT, Inc. we are making an effort to keep the planet cleaner and greener for everyone. It is our belief that we must act now to preserve our environment for future generations to come.



TABLE OF CONTENTS


| | |
|-------|--------------------------------------------------------------------|
| 04-17 | QuEChERS Line |
| 06 | QuEChERS Extraction Salts |
| 06 | QuEChERS Multi-Packs - Mylar Pouches + Centrifuge Tubes |
| 07 | Mylar Pouches - Standalone |
| 08 | Centrifuge Tube Format - Salts Pre-Loaded |
| 09 | Dispersive SPE (dSPE) Clean-up Sorbents |
| 09 | 2 ml dSPE Format |
| 10 | 15 ml dSPE Format |
| 11 | Multi-Pack - Mylar Pouches + Centrifuge Tubes |
| 11 | Mylar Pouches - Standalone |
| 12 | ChloroFiltr [®] |
| 13 | SpinFiltr [®] |
| 14-15 | LipiFiltr [®] |
| 16 | Quick QuEChERS |
| 17 | SPE Cartridge Clean-up |
| 18-23 | ENVIRO-CLEAN[®] Universal and Specialty Cartridges |
| 19 | ENVIRO-CLEAN [®] Universal Cartridges |
| 20 | ENVIRO-CLEAN [®] Specialty Cartridges |
| 21 | PFAS Solution Guide |
| 22 | ENVIRO-CLEAN [®] HL DVB |
| 23-24 | ENVIRO-CLEAN [®] 8270 |
| 25-26 | ENVIRO-CLEAN[®] Push-Through Cartridges |
| 26 | ENVIRO-CLEAN [®] Glyphosate Purification Cartridges |
| 27-29 | ENVIRO-CLEAN[®] Fusion[®] Ag+ |
| 28 | ENVIRO-CLEAN [®] Fusion [®] Ag+ Cartridges |
| 30-51 | ENVIRO-CLEAN[®] Solid-Phase Extraction Cartridges |
| 33 | Hydrophobic Phases |
| 38 | Hydrophilic Phases |
| 41 | Anion-Exchange Phases |
| 43 | Cation-Exchange Phases |
| 44 | Copolymeric Phases |
| 48 | Polymeric Phases |
| 49 | Inert Glass Cartridge |
| 50-51 | ENVIRO-CLEAN[®] Bulk Sorbent Guide |
| 51 | ENVIRO-CLEAN [®] Bulk Sorbent |
| 52-54 | HydraFlow[®] |
| 54 | HydraFlow [®] Unit Overview |
| 55-58 | Manifolds |
| 56 | 6 Station Vacuum Manifold |
| 57-58 | Glass Block Vacuum Manifold |
| 59-62 | SELECTRA[®] U/HPLC Columns |
| 60 | SELECTRA [®] DA |
| 61 | SELECTRA [®] C18 |
| 61 | SELECTRA [®] Aqueous C18 |
| 62 | SELECTRA [®] PFPP |
| 62 | SELECTRA [®] C8 |
| 63-64 | SelectraCore[®] |
| 64 | SelectraCore [®] Core-Shell Columns |
| 65-66 | SPeVAP |
| 66 | SPeVAP Specifications |
| 67-70 | Accessories |
| 68 | GC Liners |
| 69 | Reservoirs |
| 70 | Frits |
| 71 | Customer Service |

QuEChERS




QuEChERS Methods Schematic Flow Chart

STEP 1 - EXTRACTION PROCESSES



STEP 2 - DISPERSIVE SPE CLEAN-UP PROCESSES



Extraction Components and Their Functions

- MgSO₄** - Facilitates solvent partitioning
- ACN** - Provides the extraction of broadest range of pesticides with least amount of undesired matrix
- HOAc** - Adjusts pH
- Buffer Salts** - Prevent degradation of pH sensitive analytes
- PSA** - Removes sugars, acids, and some pigments
- GCB** - Strong sorbent for pigment removal
- C18** - Removes long-chain fatty compounds and other non-polar interferences

- ACN**-acetonitrile
- MgSO₄** -magnesium sulfate
- HOAc**-acetic acid
- NaCl**-sodium chloride
- Na₃Citr.**-sodium citrate tribasic dihydrate
- Na₂HCitr.**-sodium citrate dibasic sesquihydrate
- PSA**-primary secondary amine
- TPP**-triphenyl phosphate (surrogate)
- GCB**-graphitized carbon black

QuEChERS Extraction Salts

Multi-Packs - Mylar Pouches + Centrifuge Tubes

QuEChERS extraction salts for all of the popular QuEChERS methods are available in individual mylar pouches for your convenience. Each pack of 50 pouches comes with a rack of 50 empty centrifuge tubes.



| Part Number | Description | Quantity | Contents |
|-----------------|--------------------------|----------|-------------------------------------------------------------------------------------------------------------------------------------------------|
| ECMSSC50CT-MP | Non-Buffered Extraction | 50/Pack | 4000 mg MgSO ₄ 1000 mg NaCl |
| ECQUUS250CT-MP | Non-Buffered Extraction | 50/Pack | 4000 mg MgSO ₄ 2000 mg NaCl |
| ECMSSC50CTFS-MP | Non-Buffered Extraction | 50/Pack | 6000 mg MgSO ₄ 1500 mg NaCl |
| EC4MSSA50CT-MP | AOAC 2007.01 Method | 50/Pack | 4000 mg MgSO ₄ 1000 mg NaOAc |
| ECMSSA50CT-MP | AOAC 2007.01 Method | 50/Pack | 6000 mg MgSO ₄ 1500 mg NaOAc |
| ECQUEU750CT-MP | EN 15662 Method | 50/Pack | 4000 mg MgSO ₄ 1000 mg NaCl 500 mg Na ₂ HCitr.1.5H ₂ O 1000 mg Na ₃ Citr.2H ₂ O |
| EUMIV50CT-MP | EN 15662 Method | 50/Pack | 6000 mg MgSO ₄ 1500 mg NaCl 750 mg Na ₂ HCitr.1.5H ₂ O 1500 mg Na ₃ Citr.2H ₂ O |
| ECQUUS950CT-MP | THC Potency & Pesticides | 50/Pack | Proprietary Salt Blend for THC Potency & Pesticide Testing |

QuEChERS Extraction Salts

Mylar Pouch Format



| Part Number | Description | Quantity | Contents |
|-------------|-------------------------|----------|------------------------------------------------------------------------------------------------------------------------------------------------------------|
| ECMSSC-MP | Non-Buffered Extraction | 50/Pack | 4000 mg MgSO ₄ 1000 mg NaCl |
| ECQUUS2-MP | Non-Buffered Extraction | 50/Pack | 4000 mg MgSO ₄ 2000 mg NaCl |
| ECMSSCFS-MP | Non-Buffered Extraction | 50/Pack | 6000 mg MgSO ₄ 1500 mg NaCl |
| ECQUVIN-MP | Non-Buffered Extraction | 50/Pack | 8000 mg MgSO ₄ 2000 mg NaCl |
| EC4MSSA-MP | AOAC 2007.01 Method | 50/Pack | 4000 mg MgSO ₄ 1000 mg NaOAc |
| ECMSSA-MP | AOAC 2007.01 Method | 50/Pack | 6000 mg MgSO ₄ 1500 mg NaOAc |
| ECGMSSA-MP | AOAC 2007.01 Method | 50/Pack | 6000 mg MgSO ₄ (granular) 1500 mg NaOAc |
| ECQUEU7-MP | EN 15662 Method | 50/Pack | 4000 mg MgSO ₄ 1000 mg NaCl 500 mg Na ₂ HCitr.1.5H ₂ O 1000 mg Na ₃ Citr.2H ₂ O |
| ECGQUEU7-MP | EN 15662 Method | 50/Pack | 4000 mg MgSO ₄ (granular) 1000 mg NaCl 500 mg Na ₂ HCitr.1.5H ₂ O 1000 mg Na ₃ Citr.2H ₂ O |
| EUMIV-MP | EN 15662 Method | 50/Pack | 6000 mg MgSO ₄ 1500 mg NaCl 750 mg Na ₂ HCitr.1.5H ₂ O 1500 mg Na ₃ Citr.2H ₂ O |

QuEChERS Extraction Salts

Centrifuge Tube Format

UCT offers an extensive selection of QuEChERS extraction salts pre-packed in either 15 or 50 ml polypropylene centrifuge tubes.

| Part Number | Description | Volume | Quantity | Contents |
|----------------|-------------------------|--------|----------|-------------------------------------------------------------------------------------------------------------------------------------------------|
| ECMS4MSC550CT | Non-Buffered Extraction | 50 mL | 50/Pack | 4000 mg MgSO ₄ 500 mg NaCl |
| ECMSSC50CT | Non-Buffered Extraction | 50 mL | 50/Pack | 4000 mg MgSO ₄ 1000 mg NaCl |
| ECMSSC50CTFS | Non-Buffered Extraction | 50 mL | 50/Pack | 6000 mg MgSO ₄ 1500 mg NaCl |
| ECQUVIN50CT | Non-Buffered Extraction | 50 mL | 50/Pack | 8000 mg MgSO ₄ 2000 mg NaCl |
| ECMS4MSA1M50CT | AOAC 2007.01 Method | 50 mL | 50/Pack | 6000 mg MgSO ₄ 1500 mg NaOAc |
| ECMSSA50CT | AOAC 2007.01 Method | 50 mL | 50/Pack | 6000 mg MgSO ₄ 1500 mg NaOAc |
| ECQUEU750CT | EN 15662 Method | 50 mL | 50/Pack | 4000 mg MgSO ₄ 1000 mg NaCl 500 mg Na ₂ HCitr.1.5H ₂ O 1000 mg Na ₃ Citr.2H ₂ O |
| EUMIV50CT | EN 15662 Method | 50 mL | 50/Pack | 6000 mg MgSO ₄ 1500 mg NaCl 750 mg Na ₂ HCitr.1.5H ₂ O 1500 mg Na ₃ Citr.2H ₂ O |
| | | | | |
| ECQUUS1015CT | Non-Buffered Extraction | 15 mL | 50/Pack | 400 mg MgSO ₄ 100 mg NaCl |
| ECQUUS1115CT | Non-Buffered Extraction | 15 mL | 50/Pack | 800 mg MgSO ₄ 200 mg NaCl |
| ECQUUS15CT | AOAC 2007.01 Method | 15 mL | 50/Pack | 400 mg MgSO ₄ 100 mg NaOAc |
| ECQUEU415CT | EN 15662 Method | 15 mL | 50/Pack | 4000 mg MgSO ₄ 1000 mg NaCl 500 mg Na ₂ HCitr.1.5H ₂ O 1000 mg Na ₃ Citr.2H ₂ O |

Dispersive-SPE (dSPE) Sorbents

dSPE Clean-Up

UCT offers the most extensive line of dSPE products, conveniently packaged in centrifuge tube format. Each corresponding part number contains pre-weighed sorbents for use with Original, Association of Analytical Communities (AOAC) and European Committee for Standardization (CEN) QuEChERS methods. UCT dSPE kits make sample prep quick and easy for multi-residue analysis in a wide variety of matrices.

| Part Number | Volume | Quantity | Contents |
|-----------------|--------|----------|------------------------------------------------------------------|
| ECQUEU12CT | 2 mL | 100/Pack | 150 mg MgSO ₄ + 25 mg PSA |
| CUMPS2CT | 2 mL | 100/Pack | 150 mg MgSO ₄ + 50 mg PSA |
| CUMC182CT | 2 mL | 100/Pack | 150 mg MgSO ₄ + 50 mg C18 |
| ECQUEU252CT | 2 mL | 100/Pack | 150 mg MgSO ₄ + 7.5 mg GCB |
| ECQUEU22CT | 2 mL | 100/Pack | 150 mg MgSO ₄ + 25 mg PSA + 25 mg C18 |
| CUMPSC18CT | 2 mL | 100/Pack | 150 mg MgSO ₄ + 50 mg PSA + 50 mg C18 |
| CUMPS15C18CT | 2 mL | 100/Pack | 150 mg MgSO ₄ + 150 mg PSA + 50 mg C18 |
| ECQUEU32CT | 2 mL | 100/Pack | 150 mg MgSO ₄ + 25 mg PSA + 2.5 mg GCB |
| ECQUEU42CT | 2 mL | 100/Pack | 150 mg MgSO ₄ + 25 mg PSA + 7.5 mg GCB |
| CUMPSCB2CT | 2 mL | 100/Pack | 150 mg MgSO ₄ + 50 mg PSA + 50 mg GCB |
| CUMSC1875CB2CT | 2 mL | 100/Pack | 150 mg MgSO ₄ + 50 mg C18 + 7.5 mg GCB |
| ECQUUS72CT | 2 mL | 100/Pack | 50 mg PSA + 50 mg C18 + 25 mg GCB |
| CUMPSC1875CB2CT | 2 mL | 100/Pack | 150 mg MgSO ₄ + 50 mg PSA + 50 mg C18 + 7.5 mg GCB |
| ECQUEU122CT | 2 mL | 100/Pack | 150 mg MgSO ₄ + 50 mg PSA + 50 mg C18 + 50 mg GCB |

* All C18 sorbent utilized in dSPE clean-up configurations above is endcapped.

Dispersive-SPE (dSPE) Sorbents

15/50 mL Format

| Part Number | Volume | Quantity | Contents |
|-----------------|--------|----------|------------------------------------------------------------------|
| ECC18515CT | 15 mL | 50/Pack | 500 mg C18 |
| ECQUAS815CT | 15 mL | 50/Pack | 300 mg DVB |
| ECQUEU2515CT | 15 mL | 50/Pack | 450 mg MgSO ₄ + 450 mg C18 |
| ECQUAS315CT | 15 mL | 50/Pack | 750 mg MgSO ₄ + 250 mg C18 |
| ECQUUS515CT | 15 mL | 50/Pack | 900 mg MgSO ₄ + 150 mg C18 |
| ECMPS15CT | 15 mL | 50/Pack | 900 mg MgSO ₄ + 150 mg PSA |
| ECQUAS415CT | 15 mL | 50/Pack | 900 mg MgSO ₄ + 300 mg PSA |
| ECMPSA15CT | 15 mL | 50/Pack | 1200 mg MgSO ₄ + 200 mg PSA |
| ECMS12CPSA415CT | 15 mL | 50/Pack | 1200 mg MgSO ₄ + 400 mg PSA |
| ECMMCCNAX215CT | 15 mL | 50/Pack | 1200 mg MgSO ₄ + 200 mg Aminopropyl |
| ECQUEU2115CT | 15 mL | 50/Pack | 450 mg MgSO ₄ + 450 mg PSA + 150 mg C18 |
| ECQUEU2415CT | 15 mL | 50/Pack | 750 mg MgSO ₄ + 750 mg PSA + 250 mg C18 |
| ECQUEU315CT | 15 mL | 50/Pack | 900 mg MgSO ₄ + 150 mg PSA + 150 mg C18 |
| ECMPSC1815CT | 15 mL | 50/Pack | 900 mg MgSO ₄ + 300 mg PSA + 150 mg C18 |
| ECQUAS515CT | 15 mL | 50/Pack | 900 mg MgSO ₄ + 300 mg PSA + 300 mg C18 |
| CUMPSC1815CT2 | 15 mL | 50/Pack | 1200 mg MgSO ₄ + 400 mg PSA + 400 mg C18 |
| ECQUEU515CT | 15 mL | 50/Pack | 900 mg MgSO ₄ + 150 mg PSA + 15 mg GCB |
| ECQUEU615CT | 15 mL | 50/Pack | 900 mg MgSO ₄ + 150 mg PSA + 45 mg GCB |
| ECMPSCB15CT | 15 mL | 50/Pack | 900 mg MgSO ₄ + 300 mg PSA + 150 mg GCB |
| ECQUEU1115CT | 15 mL | 50/Pack | 1200 mg MgSO ₄ + 400 mg PSA + 400 mg GCB |
| ECQUHK315CT | 15 mL | 50/Pack | 700 mg MgSO ₄ + 200 mg PSA + 200 mg C18 + 35 mg GCB |
| ECQUAS615CT | 15 mL | 50/Pack | 900 mg MgSO ₄ + 450 mg PSA + 300 mg C18 + 50 mg GCB |
| ECQUUS215CT | 15 mL | 50/Pack | 1200 mg MgSO ₄ + 400 mg PSA + 400 mg C18 + 400 mg GCB |
| | | | |
| ECC1850CT | 50 mL | 50/Pack | 500 mg C18 |
| ECMSC1850CT | 50 mL | 50/Pack | 1500 mg MgSO ₄ + 500 mg C18 |

Dispersive-SPE (dSPE) Sorbents

dSPE Multi-Packs - Mylar Pouches + 15 mL Centrifuge Tubes

UCT also offers it's diverse dSPE product line in Mylar pouch format.

| Part Number | Volume | Quantity | Contents |
|--------------------|--------|----------|---------------------------------------------------|
| ECMPS15CT-MP | 15 mL | 50/Pack | 900 mg MgSO ₄ 150 mg PSA |
| ECMS12CPSA415CT-MP | 15 mL | 50/Pack | 1200 mg MgSO ₄ 400 mg PSA |
| ECMPSCB15CT-MP | 15 mL | 50/Pack | 900 mg MgSO ₄ 300 mg PSA 150 mg GCB |



dSPE Mylar Pouches

| Part Number | Quantity | Contents |
|----------------|----------|--------------------------------------------------------------|
| ECMPS-MP | 50/Pack | 900 mg MgSO ₄ + 150 mg PSA |
| ECMS12CPSA4-MP | 50/Pack | 1200 mg MgSO ₄ + 400 mg PSA |
| ECMSC18-MP | 50/Pack | 1500 mg MgSO ₄ + 500 mg C18 |
| ECQUEU315-MP | 50/Pack | 900 mg MgSO ₄ + 150 mg PSA + 150 mg C18 |
| CUMPSC182-MP | 50/Pack | 1200 mg MgSO ₄ + 400 mg PSA + 400 mg C18 |
| ECQU001-MP | 50/Pack | 750 mg MgSO ₄ + 125 mg PSA + 12.5 mg GCB |
| ECQU002-MP | 50/Pack | 750 mg MgSO ₄ + 125 mg PSA + 37.5 mg GCB |
| ECQUEU5-MP | 50/Pack | 900 mg MgSO ₄ + 150 mg PSA + 15 mg GCB |
| ECMPSCB-MP | 50/Pack | 900 mg MgSO ₄ + 300 mg PSA + 150 mg GCB |
| ECQUEU12-MP | 50/Pack | 150 mg MgSO ₄ + 50 mg PSA + 50 mg C18 + 50 mg GCB |

ChloroFiltr[®]

ChloroFiltr[®]

Traditional QuEChERS methods use graphitized carbon black (GCB) to remove chlorophyll from sample extracts. Although GCB is very effective in removing chlorophyll, it can also remove planar pesticides. UCT has developed a unique sorbent that removes chlorophyll without the loss of planar pesticides. ChloroFiltr[®] is used as an alternative to GCB during dSPE clean-up with no method modifications needed.

| Part Number | Volume | Quantity | Contents |
|---------------|--------|----------|---------------------------------------------------------------------------------------|
| CUMPSGG2CT | 2 mL | 100/Pack | 150 mg MgSO ₄ + 50 mg PSA + 50 mg Chlorofiltr [®] |
| ECQUCHL12CT | 2 mL | 100/Pack | 150 mg MgSO ₄ + 50 mg C18 + 50 mg Chlorofiltr [®] |
| CUMPSGGC182CT | 2 mL | 100/Pack | 150 mg MgSO ₄ + 50 mg PSA + 50 mg C18 + 50 mg Chlorofiltr [®] |
| ECMSGG15CT | 15 mL | 50/Pack | 900 mg MgSO ₄ + 150 mg Chlorofiltr [®] |
| ECMPSSGG15CT | 15 mL | 50/Pack | 900 mg MgSO ₄ + 300 mg PSA + 150 mg Chlorofiltr [®] |
| ECQUCHL115CT | 15 mL | 50/Pack | 900 mg MgSO ₄ + 300 mg PSA + 300 mg Chlorofiltr [®] |
| ECQUCHL215CT | 15 mL | 50/Pack | 900 mg MgSO ₄ + 300 mg C18 + 300 mg Chlorofiltr [®] |
| ECQUCHL315CT | 15 mL | 50/Pack | 900 mg MgSO ₄ + 300 mg PSA + 300 mg C18 + 300 mg Chlorofiltr [®] |
| ECMPSSGG50CT | 50 mL | 50/Pack | 1800 mg MgSO ₄ + 600 mg PSA + 300 mg Chlorofiltr [®] |
| ECQUCHL150CT | 50 mL | 50/Pack | 1800 mg MgSO ₄ + 600 mg PSA + 600 mg Chlorofiltr [®] |
| ECQUCHL250CT | 50 mL | 50/Pack | 1800 mg MgSO ₄ + 600 mg C18 + 600 mg Chlorofiltr [®] |
| ECQUCHL350CT | 50 mL | 50/Pack | 1800 mg MgSO ₄ + 600 mg PSA + 600 mg C18 + 600 mg Chlorofiltr [®] |


SpinFiltr[®]



UCT's SpinFiltr[®] takes the hassle out of dSPE

- dSPE + Ultrafiltration in a single step.
- Obtain enhanced extract purification with the built-in 0.2 μm PTFE filter
- Recover additional sample volume without worrying about the disruption of centrifugation layers.
- Simply discard dSPE chamber containing unwanted matrix and sorbet following centrifugation.

| Part Number | Top Chamber Capacity | Quantity | Contents |
|-------------|----------------------|----------|--------------------------------------------------------------------------------------|
| ECQUSF14CT | 1 mL | 50/Pack | 150 mg MgSO ₄ 50 mg PSA |
| ECQUSF24CT | 1 mL | 50/Pack | 50 mg MgSO ₄ 50 mg PSA 50 mg C18 |
| ECQUSF34CT | 1 mL | 50/Pack | 150 mg MgSO ₄ 25 mg PSA 2.5 mg GCB |
| ECQUSF44CT | 1 mL | 50/Pack | 150 mg MgSO ₄ 50 mg PSA 50 mg C18 7.5 mg GCB |
| ECQUSF54CT | 1 mL | 50/Pack | 150 mg MgSO ₄ 50 mg PSA 50 mg C18 50 mg Chlorofiltr [®] |
| ECQUSF64CT | 1 mL | 50/Pack | 150 mg MgSO ₄ 50 mg C18 50 mg Chlorofiltr [®] |
| ECQUSF74CT | 1 mL | 50/Pack | 150 mg MgSO ₄ 50 mg C18 |



Discover a whole new approach to dSPE!

Lipifiltr®

Lipifiltr™

The easy-to-use Lipifiltr® push-thru cartridge was designed to remove lipids from acetonitrile extracts. Samples are extracted using a standard QuEChERS procedure and an aliquot of the supernatant is simply pushed through the Lipifiltr® purification cartridge using a disposable syringe. The purified extract is collected in an autosampler vial and analyzed by UHPLC-MS/MS. The ability to obtain significantly cleaner extracts, ease of use, and time and cost savings make the new Lipifiltr® push-thru cartridges an attractive cleanup option for laboratories conducting pesticide residue analysis in complex fatty samples.

Lipifiltr® Push Thru Cartridge

| Part Number | Unit |
|-------------|---------|
| LPFLTR01 | 50/Pack |



LipiFiltr®

Gravimetric Analysis

The performance of the LipiFiltr® push-thru cartridges to remove lipids and other matrix components was evaluated gravimetrically by collecting 2 mL of sample before and after cleanup in pre-weighed test tubes and heating them to dryness at 110 °C in an oven.

| Matrix | Matrix Removal (%) |
|--------------|--------------------|
| Beef | 79.1 |
| Black Olives | 84.3 |
| Avocado | 54.7 |
| Salmon | 80.9 |
| Chicken fat | 71.7 |
| Olive oil | 61.5 |
| Nuts | 84.3 |
| Swordfish | 80.9 |


Extract Concentration Before and After LipiFiltr® Cleanup




Beef



Black Olives




Avocado



Salmon



Chicken Fat



Olive Oil



Nuts



Swordfish

Quick QuEChERS

Simple, Fast, Efficient Cartridges for Clean-Up of QuEChERS Extracts


UCT's QuICK QuEChERS push-thru cartridge eliminates the need for shaking and centrifugation of extracts during clean-up, significantly reducing sample processing time. In addition, any residual sorbent is filtered via the frit providing a clean, final extract for analysis.*

After QuEChERS sample extraction:

1. Draw the extract into a disposable syringe
2. Push the extract through the cartridge into a sample vial
3. Sample is ready for analysis by LC or GC

The UCT QuICK QuEChERS Cartridge provides results comparable to traditional dSPE but without the need for centrifugation.

dSPE Recovery Data vs. QuICK QuEChERS



| Part Number | Cartridge Size | Contents |
|-------------|----------------|----------------------------------------------------------------------------|
| ECPURMPSMC | Medium | 110 mg MgSO ₄ (top layer); PTFE Frit; 180 mg PSA (bottom layer) |
| ECPURMPSLC | Large | 110 mg MgSO ₄ (top layer); PTFE Frit; 600 mg PSA (bottom layer) |

* Product developed by Steven C. Moser - OK Department of Agriculture, Food and Forestry

SPE Dual Phase Cartridge Clean-Up

SPE dual phase cartridges provide an alternative clean-up option for complicated matrices and/or when dSPE does not provide adequate cleanliness of sample extracts. Examples of matrices that are suitable for SPE include tea, herbs, spices and high lipid content samples.

| Part Number | Volume | Quantity | Frit Type | Contents |
|---------------|--------|----------|-----------|-----------------------------------------------------------|
| ECPSACB6 | 6 mL | 30/Pack | PTFE | 200 mg GCB 400 mg PSA |
| ECPSACB256 | 6 mL | 30/Pack | PTFE | 250 mg GCB 500 mg PSA |
| ECPSACB256-PE | 6 mL | 30/Pack | PE | 250 mg GCB 500 mg PSA |
| ECPSACB21M6 | 6 mL | 30/Pack | PTFE | 250 mg GCB 1000 mg PSA |
| ECPSACB506 | 6 mL | 30/Pack | PTFE | 500 mg GCB 500 mg PSA |
| ECPSACB506P | 6 mL | 30/Pack | PE | 500 mg GCB 500 mg PSA |
| ECNAXCB506 | 6 mL | 30/Pack | PTFE | 500 mg GCB 500 mg CUNAX |
| EUMSPSA6 | 6 mL | 30/Pack | PTFE | 500 mg PSA 750 mg MgSO ₄ |
| ECPSAC1856 | 6 mL | 30/Pack | PTFE | 500 mg PSA 500 mg C18 |
| ECMSPSACB6 | 6 mL | 30/Pack | PTFE | 250 mg GCB 500 mg PSA 750 mg MgSO ₄ |
| EUSILMSSM26 | 6 mL | 30/Pack | PTFE | 1000 mg silica 200 mg muffled anhydrous sodium sulfate |
| | | | | |
| EUCARBC18515 | 15 mL | 20/Pack | PTFE | 500 mg GCB 500 mg C18 |
| EEMSC1811M15 | 15 mL | 20/Pack | PTFE | 1000 mg C18 1000 mg MgSO ₄ |
| EEMS2C181M15 | 15 mL | 20/Pack | PTFE | 1000 mg C18 2000 mg MgSO ₄ |
| EUPSAC181M15 | 15 mL | 20/Pack | PTFE | 1000 mg PSA 1000 mg C18 |

Enviro-Clean®


Universal & Specialty Cartridges




ENVIRO-CLEAN® Universal Cartridges




UNIVERSAL
OIL & GREASE



UNIVERSAL
525



UNIVERSAL
PAH / DRO



UNIVERSAL
C18



UNIVERSAL
549



UNIVERSAL
Zero-Blank Filter™

| Part Number | Product Name | Description | Quantity | Amount/Tube Vol. |
|-------------|------------------------------|-------------------------------------------------------------------------------|----------|------------------|
| ECUNIOGXF | UNIVERSAL OIL & GREASE | For EPA Method 1664 / Sorbent C18 | 15/Pack | 2000 mg / 83 mL |
| ECUNI525 | UNIVERSAL 525 | For EPA Method 525.2 & 525.3 / Sorbent C18 | 8/Pack | 1500 mg / 83 mL |
| ECUNIPAH | UNIVERSAL PAH / DRO | For PAH and Diesel Range Organics extractions / Sorbent C18 | 8/Pack | 2000 mg / 83 mL |
| ECUNIC18 | UNIVERSAL C18 | For extraction of pesticides, herbicides and PCBs, etc. / Sorbent C18 | 8/Pack | 1100 mg / 83 mL |
| ECUNI549 | UNIVERSAL 549 | For EPA Method 549 / Sorbent C8 | 8/Pack | 500 mg / 83 mL |
| ECUNIDVB500 | UNIVERSAL DVB | For extraction of a wide range of analytes / Sorbent PS-DVB | 8/Pack | 500 mg / 83 mL |
| ECBLANK | UNIVERSAL Zero-Blank Filter™ | Proprietary adsorbent for filtration of lab air for use during sorbent drying | 6/Pack | 200 mg / 83 mL |

ENVIRO-CLEAN® UNIVERSAL CARTRIDGE RESERVOIRS

| Part Number | Description | Quantity |
|-------------|------------------|----------|
| ERFT1FUNIP | 1 10µm PTFE frit | 10/Pack |
| ERTFT1FUNIP | 1 50µm PTFE frit | 10/Pack |
| ERFV00UNIP | Empty Reservoir | 10/Pack |

FITS THE HORIZON SPE-DEX® 4790 AUTOMATED EXTRACTION SYSTEM



ENVIRO-CLEAN[®] Specialty Cartridges

For Environmental & EPA Extractions

| Product Name | Description | Part Number | Amount/Tube Vol. | Units/Pack |
|------------------------------------------------------------------|--------------------------------------------------------------------------------------------------|--------------------------------------------|---------------------------------------------------------------|----------------|
| Clean-Elute™ | Diatomaceous Earth for EPA Method 509 | CLEAN-ELUTE | 25,000 mg / 200mL | 108 |
| Enviro-Clean® HL DVB | For use in EPA Methods 526, 528, 529, 530, 544, 553, 1694, 8321B, 8141B, 8330 | ECHLD156-P | 500 mg / 6mL | 30 |
| Enviro-Clean® 521 & 522 | Activated Carbon for EPA Methods 521 & 522 | EU52112M6 | 2000 mg / 6mL | 30 |
| Enviro-Clean® 523 | Graphitized Carbon for EPA Method 523 | EC5232506 | 250 mg / 6mL | 30 |
| Enviro-Clean® 525 | Novel C18 blend for use in EPA Method 525.2 | EC525006 | 1500 mg / 6mL | 30 |
| Enviro-Clean® 535 (90 m ² /g SA) | For EPA Method 535 or in applications requiring Graphitized Carbon Black with lower surface area | EC535156 | 500 mg / 6mL | 30 |
| Enviro-Clean® 537.1 | For EPA Method 537 Perfluorinated Alkyl Acids | ECDVB156-P | 500 mg / 6mL | 30 |
| Enviro-Clean® 541 | For EPA Method 541 Polar Organic Compounds in Drinking Water | EU541163 | 600 mg / 3mL | 50 |
| Enviro-Clean® 548 | For EPA Method 548.1 Endothall Extraction | EC548006 | 548 Slurry / 6mL | 30 |
| Enviro-Clean® C8 | For EPA Method 549 Diquat & Paraquat | EEC08156 | 500 mg / 6mL | 50 |
| Enviro-Clean® EPH Fractionation | Developed for fractionation of MA EPH (gravity flow) | XRSIHT13M15 XRSIHT15M25 CUSILHT15M25 | 3000 mg / 15mL 5000 mg / 25mL 5000 mg / 25mL | 24 20 20 |
| Enviro-Clean® Silica Gel | For silica gel clean-up applications | EUSILMSSM26 | 1000 mg silica: 200 mg muffled anhydrous sodium sulfate / 6mL | 30 |
| Enviro-Clean® Anhydrous Sodium Sulfate Drying Cartridge | Used for the removal of water from extracts prior to concentration and analysis | ECSS15M6 | 5000 mg / 6mL | 30 |
| Enviro-Clean® Anhydrous Sodium Sulfate Muffled – Glass Reservoir | Used for the removal of water from extracts prior to concentration and analysis | CUSS25M6G | 2500 mg / 6mL | 30 |
| Enviro-Clean® Alumina & Silica Dual-Phase Cartridge | For environmental clean-up | ECALNSIL25M25 | 2500 mg alumina neutral: 5000 mg silica / 25mL | 20 |
| Enviro-Clean® Sodium Sulfate & Florisil Dual-Phase Cartridge | For environmental clean-up | EUSSFL2M6 | 2000 mg sodium sulfate: 2000 mg florisil / 6mL | 30 |

PFAS Solution Guide

Enviro-Clean® WAX

UCT's polymeric weak-anion exchange (WAX) SPE cartridges feature enhanced cleanliness for minimal background in addition to unrivaled exchange capacities. Available in 100, 150, 200, and 500 mg bed sizes, these robust cartridges can be used for sample analysis in drinking water, ground water, surface water, and waste water. The use of this phase can be assessed in UCT's new US EPA Method 533 application solution in addition to an internally developed method for the sensitive quantification of 26 PFASs in drinking water, including the 14 covered in the US EPA Method 537.



Enviro-Clean® WAX
200mg and 500mg / 6mL

Enviro-Clean® HL DVB

In addition to UCT's weak-anion exchange polymeric sorbent, we also offer our ENVIRO-CLEAN® HL DVB extraction cartridge manufactured from an extremely clean, and highly cross-linked divinylbenzene based sorbent. It has been successfully used to clean up water samples for testing a wide range of analytes to include acidic, neutral (both polar & non-polar), and basic compounds. This highly retentive, reverse phase sorbent has been used successfully in US EPA Method 537 and US EPA Method 537.1.



Enviro-Clean® Free-Flowing Dual Phase Cartridges

Avoid time intensive dSPE and cumbersome dual stacked cartridges to meet the additional DOD B-15 method-required carbon purification for PFAS testing. Through the use of UCT's Enviro-Clean dual phase WAX+ Graphitized Carbon Black (GCB) cartridge, clean-up can be targeted in one step minimizing loss of long chain PFAS compounds. Cartridges are available in 2 varying configurations, 200mg of WAX+ 50mg of GCB and 500 mg of WAX+ 50mg GCB. This sample clean-up approach is quickly becoming the "go-to" method for non-potable water and soil in advance of an official EPA published method.

| EPA Method | Part Number | Description | UoM |
|-------------|---------------|---------------------------------------------------------------------------------------|--------|
| 533 | ECWAX116-P | Enviro-Clean WAX - PE Frits 100mg 6ml | 30/Pkg |
| 533 | ECWAX(150)6-P | Enviro-Clean WAX - PE Frits 150mg 6ml | 30/Pkg |
| 533 | ECWAX126-P | Enviro-Clean WAX - PE Frits 200mg 6ml | 30/Pkg |
| 533 | ECWAX156-P | Enviro-Clean WAX - PE Frits 500mg 6ml | 30/Pkg |
| 537.1 | ECHLD116-P | Enviro-Clean HL DVB - PE Frits 100mg 6ml | 30/Pkg |
| 537.1 | ECHLD(150)6-P | Enviro-Clean HL DVB - PE Frits 150mg 6ml | 30/Pkg |
| 537.1 | ECHLD126-P | Enviro-Clean HL DVB- PE Frits 200mg 6ml | 30/Pkg |
| 1633 | ECHLD156-P | Enviro-Clean HL DVB - PE Frits 500mg 6ml | 30/Pkg |
| DoD QSM 5.3 | ECWAXCB206-B | PFAS Analysis-Dual Phase Cartridge-200mg ECWAX+50mg Graphitized Carbon Black-PE Frits | 30/Pkg |
| DoD QSM 5.3 | ECWAXCB506-B | PFAS Analysis-Dual Phase Cartridge-500mg ECWAX+50mg Graphitized Carbon Black-PE Frits | 30/Pkg |


ENVIRO-CLEAN[®] Specialty Cartridges

ENVIRO-CLEAN[®] HL DVB

ENVIRO-CLEAN[®] HL DVB extraction columns are manufactured from an extremely clean, highly cross-linked divinylbenzene based sorbent. The material was developed with the environmental market in mind. It has been successfully used to extract a wide range of analytes from water samples. By varying the sample pH, wash and elution solvents ENVIRO-CLEAN[®] HL DVB can be used to analyze acidic, basic, and neutral (both polar & non-polar) compounds.



| Part Number | Description | Quantity |
|---------------|------------------------|----------|
| ECHLD(150)6-P | 150 mg / 6mL Cartridge | 30/Pack |
| ECHLD126-P | 200 mg / 6mL Cartridge | 30/Pack |
| ECHLD156-P | 500 mg / 6mL Cartridge | 30/Pack |



ENVIRO-CLEAN[®] Specialty Cartridges

EPA METHOD 8270 - Extraction of Acids, Bases, and Neutrals in Water using Solid Phase Extraction

UCT offers a 2 cartridge system and extraction procedure for EPA Method 8270. A wide range of 137 target analytes and 6 surrogates can be successfully analyzed using this method. The procedure is reliable, efficient, and cost-effective. The tandem cartridge system uses UCT's proprietary 8270 cartridge in-line with our activated carbon cartridge. High throughput can be achieved by extracting multiple samples simultaneously using a multi-port SPE manifold combined with a 12 position collection rack, which allows for the simultaneous extraction of up to 12 samples at once. A set of 24 samples can be extracted in 5 to 6 hours.

Product Benefits

- Cost-effective
- Reduced usage of organic solvents
- Only one sample pass is needed

5-6 hrs for a batch of 24 samples.

- No emulsion or white precipitate generated
- Shorter solvent evaporation time
- Shorter sample turnaround time
- High sample throughput
- Excellent recovery
- Cleaner extracts and chromatograms
- Cartridge body manufactured from special, proprietary polypropylene – minimizing potential source of interferences
- Packaged in Mylar to maintain cleanliness



ENVIRO-CLEAN[®] Specialty Cartridges

8270 Cartridge Kits

| 1 Liter Sample Size | | |
|---------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------|
| Part Number | Description | Quantity |
| EC8270-KIT1L | ENVIRO-CLEAN [®] 8270 STARTER KIT | Kit |
| Contents | 30 x 8270 Extraction Cartridges (p/n EC82702M15), 30 x Carbon Extraction Cartridges (p/n EU52113M6), 30 x Cartridge Adapters (p/n AD0000AS), 12 x Large Volume Transfer Tubes (p/n VMFSTFR12) | |
| EC8270-1000REFL | ENVIRO-CLEAN [®] 8270 REFILL KIT | Kit |
| Contents | 30 x 8270 Extraction Cartridges (p/n EC82702M15), 30 x Carbon Extraction Cartridges (p/n EU52113M6) | |
| 500 mL Sample Size | | |
| Part Number | Description | Quantity |
| EC8270-KIT | ENVIRO-CLEAN [®] 8270 STARTER KIT | Kit |
| Contents | 30 x 8270 Extraction Cartridges (p/n EC82701M15), 30 x Carbon Extraction Cartridges (p/n EU52112M6), 30 x Cartridge Adapters (p/n AD0000AS), 12 x Large Volume Transfer Tubes (p/n VMFSTFR12) | |
| EC8270-500REFL | ENVIRO-CLEAN [®] 8270 REFILL KIT | Kit |
| Contents | 30 x 8270 Extraction Cartridges (p/n EC82701M15), 30 x Carbon Extraction Cartridges (p/n EU52112M6) | |
| <100 mL Sample Size | | |
| Part Number | Description | Quantity |
| EC8270-KIT | ENVIRO-CLEAN [®] 8270 STARTER KIT | Kit |
| Contents | 30 x 8270 Extraction Cartridges (p/n EC8270506), 30 x Carbon Extraction Cartridges (p/n EU5211M6), 30 x Cartridge Adapters (p/n AD0000AS), 12 x Large Volume Transfer Tubes (p/n VMFSTFR12) | |
| EC8270-500REFL | ENVIRO-CLEAN [®] 8270 REFILL KIT | Kit |

| 8270 Cartridge Refills | | |
|------------------------|-----------------------------|----------|
| Part Number | Description | Quantity |
| EC8270506 | 500 mg in a 6mL Cartridge | 30/Pack |
| EC82701M15 | 1000 mg in a 15mL Cartridge | 30/Pack |
| EC82702M15 | 2000 mg in a 15mL Cartridge | 30/Pack |

| Carbon Cartridge Refills | | |
|--------------------------|----------------------------|----------|
| Part Number | Description | Quantity |
| EU5211M6 | 1000 mg in a 6mL Cartridge | 30/Pack |
| EU52112M6 | 2000 mg in a 6mL Cartridge | 30/Pack |
| EU52113M6 | 3000 mg in a 6mL Cartridge | 30/Pack |

| Extraction Manifold | | |
|---------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------|
| Part Number | Description | Quantity |
| VMF016GL | 16 Position Complete Vacuum Manifold System | Complete Unit |
| Contents | 1 x Glass Block, 1 x 16 Position Corian Lid 1 x Cover Gasket, 1 x Vacuum Gauge 1 x 16 Position Adjustable Collection Rack 1 x Glass Block Safety Tray, 16 x PTFE Tips 16 x Bulkhead Luer Fittings, 16 x Plugs | |

ENVIRO-CLEAN® Push-Thru Cartridges

ENVIRO-CLEAN® Push-Thru Purification Cartridges

UCT's wide sorbent range is also offered in a convenient push-thru format, providing simple, fast, and efficient sample clean-up. The analyst can choose from a wide range of sorbent types tailored to the specific requirements of each analysis. The small, medium, and large push-thru cartridges allow for the filtration and removal of unwanted matrix that can otherwise lead to significant ion suppression during sample analysis.



SMALL



MEDIUM




LARGE

| Part Number | Cartridge Size | Quantity | Frit Type | Contents |
|-------------|----------------|----------|-----------|-------------------------------------------------|
| EEC18MC | Medium | 50/Pack | PTFE | Enviro-Clean® Encapped C18 |
| EEC1815MC-P | Medium | 50/Pack | PE | Enviro-Clean® Encapped C18 |
| EUCARBMC | Medium | 50/Pack | PTFE | Enviro-Clean® Graphitized Carbon – 120/400 Mesh |
| EUFLSLC | Large | 50/Pack | PTFE | Enviro-Clean® Florisil® PR |
| EUALNLC-P | Large | 50/Pack | PE | Enviro-Clean® Alumina Neutral |
| EUQAXLC | Large | 50/Pack | PTFE | Enviro-Clean® Quaternary Amine |
| EUPSAMC | Medium | 50/Pack | PTFE | Enviro-Clean® Primary/Secondary Amine |
| EUPSALC | Large | 50/Pack | PTFE | Enviro-Clean® Primary/Secondary Amine |

ENVIRO-CLEAN® Push-Thru Cartridges

ENVIRO-CLEAN® Push-Thru Glyphosate Purification Cartridges

Remove unwanted matrix interferences that can lead to significant suppression and loss of ionization for glyphosate and glufosinate. Water and vegetation/fruits/swabs and other non-soil samples extracted with water are simply pushed through the Glyphosate cartridge to purify final extracts.



| Part Number | Size |
|---------------------------------------|-----------------|
| ECGLYSC | Small Cartridge |
| ECGLYLC | Large Cartridge |
| * Proprietary polymeric sorbent blend | |

* Source: Steven C. Moser OK Department of Agriculture, Food & Forestry

**Glyphosate is found at 168 m/z in ESI- mode.


Infused samples were injected onto LC/MS Ion Trap in full scan mode for an average of a 30 second infusion pre and post clean-up using UCT's Glyphosate Purification Cartridges. Background matrix peaks that can lead to significant suppression and compete with glyphosate and glufosinate during ionization were significantly reduced following clean-up.

ENVIRO-CLEAN[®] Fusion[®] Ag⁺

Fractionation of aliphatic and aromatic hydrocarbons are routinely assessed to identify the risks posed by petroleum in the environment. Solid Phase Extraction (SPE) with heat-treated silica gel is the traditional approach for the fractionation of aliphatic and aromatic hydrocarbons.

The biggest challenge of this approach is the deactivation of the silica gel due to its hygroscopic nature, which leads to inconsistent results, incomplete separation, and low recoveries. The process to optimize each batch of silica cartridges is time consuming and exhausts copious amounts of hexane solvent.

UCT's Enviro-Clean Fusion[®] Ag⁺ SPE sorbent is designed to help overcome the challenges associated with traditional silica gel fractionation. The new EPH fractionation sorbent consists of silver ions anchored onto a solid support to provide the best selectivity, capacity, and performance.



Aromatic hydrocarbons form a charge-transfer complex with silver ions ensuring complete separation of aliphatic from aromatic fractions without breakthrough.






Figure 2. Illustration of Ag⁺ interaction with aromatic and aliphatic hydrocarbons.

ENVIRO-CLEAN[®] Fusion[®] Ag⁺

ENVIRO-CLEAN[®] Fusion[®] Ag⁺ Order Information

| Description | Part Number |
|------------------------------------------------------------|-------------|
| Fusion [®] Ag ⁺ / 0.5 g, 6mL cartridge | ECFUSAG156 |
| Fusion [®] Ag ⁺ / 1g, 6mL cartridge | ECFUSAG1M6 |
| Fusion [®] Ag ⁺ / 2g, 6mL cartridge | ECFUSAG2M6 |





Experiment conditions: Fusion[®] Ag⁺ 1g/6mL cartridge; Sample: NJDEP EPH 10/08 Rev.2 aliphatics calibration standard and aromatic calibration standard; Concentration: 200 µg/mL in hexane; GC-MS system: Thermo Trace 1300 GC & ISQ MS; Restek Rxi[®]-5sil MS, 30m x 0.25mm, 0.25µm with Integra-Guard; Injection volume: 1µL split (1:100) at 300°C; GC liner: 4 mm split liner with deactivated glass wool; Temperature: Transfer line = 275°C; Ion source = 275°C; Oven temp. program: Initial temperature at 50°C, hold for 3 min; ramp at 10°C/min to 320°C, hold for 10 min; Full scan range: 35-600 amu.

ENVIRO-CLEAN[®] Fusion[®] Ag⁺

Fusion[®] Ag⁺ requires less solvent consumption and reduced extraction times as compared to silica cartridges with a similar configuration (Table 2, 3). Due to the strong ionic interaction of the silver ions with the aromatic hydrocarbon fraction, there is no need to optimize the hexane volume for each lot of silica cartridges. This results in a simplified process procedure.

Fewer solvent requirements and shortened extraction time improves cost efficiencies and laboratory productivity. The small configurations: 0.5g/6mL and 1g/6mL are ideal for automated high throughput analysis. The 2g/6mL configuration is best for higher loading capacity. In addition, using acetone instead of dichloromethane as the elution solvent for the aromatic fraction is more environmentally friendly.



Experimental conditions are the same as Figure 3. Three samples were tested from each lot.

Table 2 - Fusion[®] Ag⁺ reduces solvent consumption

| Solvent Consumption (mL) | | | |
|--------------------------------------------|-----------------|--------------------------------------------|----------------------------|
| Extraction Steps | | Fusion [®] Ag ⁺ 1g/6mL | Heat-treated Silica 2g/6mL |
| Conditioning | Acetone | 4 | |
| | Dichloromethane | | 3 |
| | Hexane | 4 | 6 |
| Sample Loading | Hexane | 1 | 1 |
| Aliphatic Fraction | Hexane | 4 | 4 |
| | Acetone | 4 | |
| Aromatic Fraction | Acetone | 4 | |
| | Dichloromethane | | 6 |
| Total Solvent Consumption (mL) | | 9 | 20 |
| Total Chlorinated Solvent Used (mL) | | 0 | 9 |

Table 3 - Fusion[®] Ag⁺ shortens extraction time

| Extraction Time (min) | | | |
|------------------------------------|------------------------------|--------------------------------------------|----------------------------|
| Extraction Steps | | Fusion [®] Ag ⁺ 1g/6mL | Heat-treated Silica 2g/6mL |
| Conditioning | Acetone | 1.8 | |
| | Dichloromethane (5-min hold) | | 5.8 |
| | Hexane | 1.7 | 1.9 |
| Sample Loading | Hexane | 0.5 | 0.5 |
| Aliphatic Fraction | Hexane | 1.7 | 1.4 |
| | Acetone | 1.5 | |
| Aromatic Fraction | Acetone | 1.5 | |
| | Dichloromethane | | 2.8 |
| Total Extraction Time (min) | | 7.2 | 12.4 |

Enviro-Clean®

Solid-Phase Extraction Cartridges



**HYDROPHOBIC / HYDROPHILIC /
ION EXCHANGE / COPOLYMERIC / POLYMERIC**

ENVIRO-CLEAN[®] Sorbents

ENVIRO-CLEAN[®] solid-phase extraction (SPE) cartridges are designed specifically for the isolation and purification of environmental analytes such as pesticides, herbicides, polyaromatic hydrocarbons, polychlorinated biphenyls and other environmentally related compounds. By utilizing ultra-clean extraction sorbents along with chemically resistant PTFE frits, an end-user not only has the ability to purify complex matrices, but also reduce ion suppression or enhancement, and most importantly enrich compounds present at trace concentration levels.

To successfully conduct SPE, a mechanistic understanding of the interaction between sorbent and analyte of interest is vital for producing optimal results. The most common retention mechanisms include non-polar interactions (van der Waals forces), polar interactions (hydrogen bonding, dipole-dipole forces), and ionic interactions (cation-anion exchange).

Non-polar phases are universal but often are considered to exhibit the least selective retention mechanisms when compared to normal phase or ion-exchange SPE. C18 is the most widely used of these phases. While this retention mechanism happens to lack specificity, it is very useful for extracting analytes that are very diverse in structure within the same sample. Several EPA approved methods for analyzing organics in water mandate the use of a C18 phase.

Retention of an analyte under normal phase conditions is primarily due to interactions between polar functional groups of the analyte and polar groups on the sorbent surface. This could include hydrogen bonding, pi-pi, and dipole-dipole interactions. This mode

is classically used to separate neutral organic compounds whose chemical nature ranges from hydrophobic to moderately polar.

Ion-exchange phases are often applied when analytes of interest carry a charge while in solution. The primary retention mechanism of the compound is based mainly on the electrostatic attraction of the charged functional group on the compound to the charged group that is bonded to the sorbent surface. ENVIRO-CLEAN[®] sorbents are available in either cation or anion exchangers and exhibit both weak and strong characteristics.

Lastly, copolymeric phases offer a dual retention mechanism providing superior cleanliness. Hydrophobic interactions in addition to ion exchange contribute to a higher degree of analyte selectivity than what was previously possible. Compounds of interest can be retained by multiple mechanisms, resulting in greater removal of matrix-related contamination.

SPE TERMINOLOGY

- Sorbent – The solid-phase material to which analytes attach during the extraction process
- Bed Volume – The quantity of solvent needed to cover the sorbent
- Capacity – The amount of analyte that a sorbent can retain
- Activation – a process of rinsing the sorbent with a solvent to clean the bed and extend the bonded groups (i.e. C18) to maximize its effectiveness
- Wash Solvent – Solvent used to wash interferences off of the sorbent prior to elution
- Retention – The attraction a solid-phase has for the analyte that causes the analyte to “adsorb” to the sorbent
- Elution – The process of removing an analyte from the sorbent for analysis

ENVIRO-CLEAN[®] SPE Cartridges

HOW TO READ ENVIRO-CLEAN[®] PART NUMBERS: ENVIRO-CLEAN[®] CODE

EU Unencapped Cartridge
EE: Encapped Cartridge

BCX Sorbent Type

153 Cartridge Size (mL)
Packing Size (x 100mg)
("M" signifies x 1,000mg)


SINGLE & DUAL PHASE
1: Cartridge contains given sorbent type.
2: Cartridge contains given sorbent type plus a hydrophobic phase.

EUBCX153:
ENVIRO-CLEAN[®] Unencapped Benzenesulfonic Acid Cation Exchange Cartridge, Single Phase, 500mg/3mL.

| Sorbent Type | |
|---------------|---------------------------------------|
| Code | Description |
| C08, C18, C30 | Carbon Chains |
| SIL | Unbonded Silica |
| PSA | n-2 Aminoethyl |
| BCX | Benzenesulfonic Acid Cation Exchanger |
| PCX | Propylsulfonic Acid Cation Exchanger |
| CCX | Carboxylic Acid Cation Exchanger |
| QAX | Quaternary Amine Anion Exchanger |
| NAX | Aminopropyl Anion Exchanger |
| FLS | Florisil [®] PR |
| ALA | Alumina - Acid |
| ALB | Alumina - Base |
| ALN | Alumina - Neutral |
| CNP | Cyanopropyl |
| CYH | Cyclohexyl |
| DOL | Diol |
| PHY | Phenyl |

HYDROPHOBIC EXTRACTION SORBENTS


UCT's ENVIRO-CLEAN[®] hydrophobic sorbents feature silanol groups at the surface of the raw silica packing that have been chemically modified with hydrophobic alkyl or aryl functional groups. These phases are commonly utilized to extract compounds that exhibit medium non-polar characteristics from a variety of complex matrices. The C18 phase is the most widely used phase for non-polar interactions because of its non-selective nature; C18 will extract a large number of compounds with differing chemical properties. To enhance selectivity, UCT offers a variety of hydrophobic sorbents ranging from C2 all the way to C30, as well as endcapped and unendcapped versions.





ENVIRO-CLEAN[®] Hydrophobic Phase

MECHANISM OF HYDROPHOBIC BONDING

Retention of organic analytes from polar solutions (e.g. water) onto these SPE materials is due primarily to the attractive forces between the carbon-hydrogen bonds in the analyte and the functional groups. These nonpolar-nonpolar attractive forces are commonly called van der Waals forces, or dispersion forces. To elute an adsorbed compound from a reversed phase sorbent, use a nonpolar solvent to disrupt the forces that bind the compound to the packing. Some polar solvents, such as methanol and acetonitrile have enough non-polar characteristics to disrupt nonpolar binding triggering compound elution. Methanol can be used as well, although it should be noted that it will remove both polar and non-polar analytes of interest as well as interferences.





Hydrophobic Sorbents & Structures



| Sorbent | Structure |
|---------------|--------------------------------------------------------------------------------------------|
| C2 Ethyl | -SiCH ₂ CH ₃ |
| C4 n-Butyl | -Si(CH ₂) ₃ CH ₃ |
| C8 Octyl | -Si(CH ₂) ₇ CH ₃ |
| C18 Octadecyl | -Si(CH ₂) ₁₇ CH ₃ |
| C30 Tricontyl | -Si(CH ₂) ₂₉ CH ₃ |
| Cyclohexyl | -Si-  |
| Phenyl | -Si-  |

ENDCAPPED VS. UNENDCAPPED

Bonded phases are manufactured by the reaction of organosilanes with activated silica. During the polymerization reaction of carbon chains to the silica backbone, a very stable silyl ether linkage forms. Our unendcapped columns allow hydroxyl sites to remain, thus making these columns slightly hydrophilic. In order to minimize this slight polarity, the hydroxyl sites are deactivated. Proprietary bonding techniques ensure that these sites are 100% reacted, leading to complete endcapping. Because there are no hydroxyl sites left, our endcapped columns render the silica non-acidic and non-polar.



ENVIRO-CLEAN[®] Hydrophobic Phase



C18, OCTADECYL

Organic Loading = 21.5%
Surface Area = 500 m²/g

Average Pore Size = 60Å
Pore Volume = 0.77 cm³/g

C8, OCTYL

Organic Loading = 11.1%
Surface Area = 500 m²/g

Average Pore Size = 60Å
Pore Volume = 0.77 cm³/g

| COLUMNS | | | | |
|------------------|---------------------|----------------|-----------|-------------|
| Tube Volume (mL) | Sorbent Amount (mg) | Units per Pack | Endcapped | Part Number |
| 1 | 100 | 100 | NO | EUC18111 |
| 1 | 100 | 100 | YES | EEC18111 |
| 3 | 200 | 50 | NO | EUC18123 |
| 3 | 200 | 50 | YES | EEC18123 |
| 3 | 500 | 50 | NO | EUC18153 |
| 3 | 500 | 50 | YES | EEC18153 |
| 6 | 500 | 50 | NO | EUC18156 |
| 6 | 500 | 50 | YES | EEC18156 |
| 6 | 1000 | 30 | NO | EUC181M6 |
| 6 | 1000 | 30 | YES | EEC181M6 |
| 10 | 100 | 50 | YES | EEC1811Z |
| 10 | 200 | 50 | NO | EUC1812Z |
| 10 | 200 | 50 | YES | EEC1812Z |
| 10 | 500 | 50 | NO | EUC1815Z |
| 10 | 500 | 50 | YES | EEC1815Z |
| 15 | 1000 | 20 | YES | EEC1811M15 |
| 15 | 2000 | 20 | NO | EUC1812M15 |
| 15 | 2000 | 20 | YES | EEC1812M15 |
| 25 | 5000 | 20 | NO | EUC1815M25 |
| 25 | 5000 | 20 | YES | EEC1815M25 |
| 75 | 10000 | 10 | NO | EUC18110M75 |

| COLUMNS | | | | |
|------------------|---------------------|----------------|-----------|-------------|
| Tube Volume (mL) | Sorbent Amount (mg) | Units per Pack | Endcapped | Part Number |
| 1 | 50 | 100 | NO | EUC081L1 |
| 1 | 100 | 100 | NO | EUC08111 |
| 1 | 100 | 100 | YES | EEC08111 |
| 3 | 200 | 50 | NO | EUC08123 |
| 3 | 200 | 50 | YES | EEC08123 |
| 3 | 500 | 50 | NO | EUC08153 |
| 3 | 500 | 50 | YES | EEC08153 |
| 6 | 500 | 50 | NO | EUC08156 |
| 6 | 500 | 50 | YES | EEC08156 |
| 6 | 1000 | 30 | NO | EUC081M6 |
| 6 | 1000 | 30 | YES | EEC081M6 |
| 10 | 100 | 50 | NO | EUC0811Z |
| 10 | 200 | 50 | NO | EUC0812Z |
| 10 | 200 | 50 | YES | EEC0812Z |
| 10 | 500 | 50 | NO | EUC0815Z |
| 10 | 500 | 50 | YES | EEC0815Z |
| 15 | 2000 | 20 | NO | EUC0812M15 |
| 15 | 2000 | 20 | YES | EEC0812M15 |
| 25 | 5000 | 20 | NO | EUC0815M25 |
| 25 | 2000 | 20 | YES | EEC0815M25 |
| 75 | 5000 | 20 | NO | EUC08110M75 |

ENVIRO-CLEAN[®] Hydrophobic Phase

C2, ETHYL

Organic Loading = 6.2%
Surface Area = 500 m²/g

Average Pore Size = 60Å
Pore Volume = 0.77 cm³/g

| COLUMNS | | | | |
|------------------|---------------------|----------------|-----------|-------------|
| Tube Volume (mL) | Sorbent Amount (mg) | Units per Pack | Endcapped | Part Number |
| 6 | 500 | 50 | YES | EEC02156 |
| 10 | 500 | 50 | YES | EEC0215Z |

C30, TRICONTYL

Organic Loading = 20.0%
Surface Area = 500 m²/g

Average Pore Size = 60Å
Pore Volume = 0.77 cm³/g

| COLUMNS | | | | |
|------------------|---------------------|----------------|-----------|-------------|
| Tube Volume (mL) | Sorbent Amount (mg) | Units per Pack | Endcapped | Part Number |
| 3 | 200 | 50 | NO | EUC30123 |
| 3 | 500 | 50 | NO | EUC30153 |
| 6 | 1000 | 30 | YES | EEC301M6 |
| 10 | 500 | 50 | NO | EUC3015Z |

PHY, PHENYL

Organic Loading = 10.8%
Surface Area = 500 m²/g

Average Pore Size = 60Å
Pore Volume = 0.77 cm³/g

| COLUMNS | | | | |
|------------------|---------------------|----------------|-----------|-------------|
| Tube Volume (mL) | Sorbent Amount (mg) | Units per Pack | Endcapped | Part Number |
| 3 | 500 | 50 | NO | EUPHY153 |
| 3 | 500 | 50 | YES | EEPHY153 |

CYH, CYCLOHEXYL


Organic Loading = 11.6%
Surface Area = 500 m²/g

Average Pore Size = 60Å
Pore Volume = 0.77 cm³/g

| COLUMNS | | | | |
|------------------|---------------------|----------------|-----------|-------------|
| Tube Volume (mL) | Sorbent Amount (mg) | Units per Pack | Endcapped | Part Number |
| 3 | 500 | 50 | NO | EUCYH153 |
| 6 | 500 | 50 | YES | EECYH156 |

ENVIRO-CLEAN[®] HYDROPHILIC NORMAL PHASE EXTRACTION SORBENTS

Hydrophilic sorbents are composed of a silica backbone bonded with carbon chains containing polar functional groups. Examples of phases that feature this functionality include bare silica, diol, and cyanopropyl phases.



ENVIRO-CLEAN[®] Hydrophobic Phase


MECHANISM OF HYDROPHILIC BONDING

Compounds are retained on hydrophilic sorbents through polar interactions including hydrogen bonding, pi-pi or dipole-dipole interactions. These types of interactions occur when the distribution of electrons between individual atoms in functional groups is unequal, causing negative and positive polarity. Polar-functionalized bonded silicas and polar adsorption media are typically used under normal phase conditions.

Hydrophilic Sorbents & Structures

| Sorbent | Structure |
|-------------|--------------------------------------------------------------------------|
| Silica | -SiOH |
| Diol | -Si(CH ₂) ₃ OCH ₂ OHCH ₂ OH |
| Cyanopropyl | -Si(CH ₂) ₃ CN |

Example of Hydrophilic Bonding



UNBONDED SILICA, ACID WASHED

Organic Loading = N/A
Surface Area = 500 m²/g

Average Pore Size = 60Å
Pore Volume = 0.77 cm³/g

COLUMNS

| Tube Volume (mL) | Sorbent Amount (mg) | Units per Pack | Part Number |
|------------------|---------------------|----------------|-------------|
| 1 | 100 | 100 | EUSIL111 |
| 3 | 200 | 50 | EUSIL123 |
| 3 | 500 | 50 | EUSIL153 |
| 6 | 500 | 50 | EUSIL156 |
| 6 | 1000 | 30 | EUSIL1M6 |
| 15 | 2000 | 20 | EUSIL12M15 |
| 25 | 5000 | 20 | EUSIL15M25 |

PHARMA-SIL[®]

Organic Loading = N/A
Surface Area = 500 m²/g

Average Pore Size = 60Å
Pore Volume = 0.82 cm³/g

COLUMNS

| Tube Volume (mL) | Sorbent Amount (mg) | Units per Pack | Part Number |
|------------------|---------------------|----------------|-------------|
| 1 | 100 | 100 | EPHSIL111 |
| 3 | 500 | 50 | EPHSIL153 |
| 6 | 500 | 50 | EPHSIL156 |
| 6 | 1000 | 30 | EPHSIL1M6 |
| 10 | 200 | 50 | EPHSIL12Z |
| 15 | 2000 | 20 | EPHSIL12M15 |
| 25 | 5000 | 20 | EPHSIL15M25 |

ALUMINA, ACIDIC

COLUMNS

| Tube Volume (mL) | Sorbent Amount (mg) | Units per Pack | Part Number |
|------------------|---------------------|----------------|-------------|
| 3 | 200 | 50 | EUALA123 |
| 3 | 500 | 50 | EUALA153 |
| 6 | 500 | 50 | EUALA156 |
| 6 | 1000 | 30 | EUALA1M6 |
| 6 | 2000 | 30 | EUALA2M6 |

ALUMINA, BASIC

COLUMNS

| Tube Volume (mL) | Sorbent Amount (mg) | Units per Pack | Part Number |
|------------------|---------------------|----------------|-------------|
| 1 | 50 | 100 | EUALB1L1 |
| 3 | 500 | 50 | EUALB153 |
| 6 | 500 | 50 | EUALB156 |
| 6 | 1000 | 30 | EUALB1M6 |

ENVIRO-CLEAN[®] Hydrophilic Phase

ALUMINA, NEUTRAL

| COLUMNS | | | |
|------------------|---------------------|----------------|-------------|
| Tube Volume (mL) | Sorbent Amount (mg) | Units per Pack | Part Number |
| 1 | 50 | 100 | EUALN1L1 |
| 1 | 100 | 100 | EUALN111 |
| 3 | 200 | 50 | EUALN123 |
| 6 | 500 | 50 | EUALN156 |
| 6 | 1000 | 30 | EUALN1M6 |
| 15 | 2000 | 20 | EUALN12M15 |
| 25 | 5000 | 20 | EUALN15M25 |

CARBON, GRAPHITIZED NON-POROUS, 120/400 MESH

| COLUMNS | | | |
|------------------|---------------------|----------------|-------------|
| Tube Volume (mL) | Sorbent Amount (mg) | Units per Pack | Part Number |
| 1 | 100 | 100 | EUCARB111 |
| 3 | 200 | 50 | EUCARB123 |
| 3 | 500 | 50 | EUCARB153 |
| 6 | 200 | 50 | EUCARB126 |
| 6 | 250 | 30 | EUCARB2L6 |
| 6 | 500 | 50 | EUCARB156 |
| 6 | 1000 | 30 | EUCARB1M6 |
| 10 | 100 | 50 | EUCARB11Z |
| 10 | 200 | 50 | EUCARB12Z |



DIOL

Organic Loading = 8.0%
Surface Area = 500 m²/g

Average Pore Size = 60 Å
Pore Volume = 0.77 cm³/g

| COLUMNS | | | |
|------------------|---------------------|----------------|-------------|
| Tube Volume (mL) | Sorbent Amount (mg) | Units per Pack | Part Number |
| 1 | 50 | 100 | EUDOL1L1 |
| 1 | 100 | 100 | EUDOL111 |
| 1 | 200 | 100 | EUDOL121 |
| 3 | 200 | 50 | EUDOL123 |
| 3 | 500 | 50 | EUDOL153 |
| 6 | 500 | 50 | EUDOL156 |
| 6 | 1000 | 30 | EUDOL1M6 |
| 10 | 500 | 50 | EUDOL15Z |
| 15 | 2000 | 20 | EUDOL12M15 |

FLORISIL[®] PR

Florisil[®] is the trademark of U.S. Silica Co.

| COLUMNS | | | |
|------------------|---------------------|----------------|-------------|
| Tube Volume (mL) | Sorbent Amount (mg) | Units per Pack | Part Number |
| 1 | 100 | 100 | EUFLS111 |
| 3 | 200 | 50 | EUFLS123 |
| 3 | 500 | 50 | EUFLS153 |
| 6 | 500 | 50 | EUFLS156 |
| 10 | 100 | 50 | EUFLS11Z |
| 10 | 500 | 50 | EUFLS15Z |
| 15 | 2000 | 20 | EUFLS12M15 |
| 25 | 5000 | 20 | EUFLS15M25 |

FLORISIL[®] A (100 / 120 Mesh)

Florisil[®] is the trademark of U.S. Silica Co.



| COLUMNS | | | |
|------------------|---------------------|----------------|-------------|
| Tube Volume (mL) | Sorbent Amount (mg) | Units per Pack | Part Number |
| 3 | 500 | 50 | EUFLSA153 |

ENVIRO-CLEAN[®] Ion Exchange Phase

MECHANISM OF ION EXCHANGE BONDING

The retention mechanism in ion exchange bonding is the electrostatic attraction of the charged functional group on the compound to the charged functional group on the SPE sorbent. In order for an optimal interaction, both the compound of interest and the functional group on the bonded silica must be fully charged. To ensure 99% or more ionization, the pH should be at least two pH units below the pKa of the cation and two pH units above the pKa of the anion. Elution occurs by using a solution/buffer to raise the pH above the pKa of the cationic group or to lower the pH below the pKa of the anion to disrupt retention. At this point, the sorbent or compound is neutralized. Ion exchange resins come in two varying forms: strong and weak. The number of charges on a strong ion exchanger remains constant regardless of the buffer pH. Weak ion exchangers display pH-dependent functionality and therefore deliver optimal performance over only a small pH range.

| Functionality | Ionization | pH units away from pKa | | | | |
|---------------|--------------|------------------------|---------|--------|---------|---------|
| | | 2 < pKa | 1 < pKa | At pKa | 1 > pKa | 2 > pKa |
| ACID | Anionic (-) | 1 | 9 | 50 | 91 | 99 |
| BASE | Cationic (+) | 99 | 91 | 50 | 9 | 1 |





ENVIRO-CLEAN[®] Ion Exchange Phase

ION EXCHANGE SORBENTS & STRUCTURES

| Sorbent | Structure | pKa |
|----------------------------------------|---------------------------------------------------------------------------------------------------------------------|----------------|
| Anion Exchangers | | |
| NAX - Aminopropyl (1° amine) | -Si-(CH ₂) ₃ NH ₂ | 9.8 |
| PSA - N-2 Aminoethyl (1° & 2° amine) | -Si-(CH ₂) ₃ NH(CH ₂) ₂ NH ₂ | 10.1, 10.9 |
| DAX - Diethylamino (3° amine) | -Si-(CH ₂) ₃ N(CH ₂ CH ₃) ₂ | 10.6 |
| EUQAX - Quaternary Amine Chloride | -Si-(CH ₂) ₃ N ⁺ (CH ₃) ₃ Cl ⁻ | Always charged |
| EHQAX - Quaternary Amine Hydroxide | -Si-(CH ₂) ₃ N ⁺ (CH ₃) ₃ OH ⁻ | Always charged |
| EAQAX - Quaternary Amine Acetate | -Si-(CH ₂) ₃ N ⁺ (CH ₃) ₃ CH ₃ COO ⁻ | Always charged |
| EFQAX - Quaternary Amine Formate | -Si-(CH ₂) ₃ N ⁺ (CH ₃) ₃ HCOO ⁻ | Always charged |
| PAX - Polyimine | -Si-(CH ₂) ₃ -R [NHCH ₃ CH ₃] _x | |

Cation Exchangers

| | | |
|---------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------|----------------|
| CCX - Carboxylic Acid | -Si-CH ₂ COOH | 4.8 |
| PCX - Propylsulfonic Acid | -Si-(CH ₂) ₃ SO ₃ H | <1 |
| BCX - Benzenesulfonic Acid | -Si-(CH ₂) ₂  ⁻ | Always charged |
| BCXHL- Benzenesulfonic Acid High Load | -Si-(CH ₂) ₂  ⁻ | Always charged |
| TAX - Triacetic Acid | -Si-(CH ₂) ₃ NH-(CH ₂) ₂ N(CH ₂ COOH) ₂ CH ₂ COOH | |

| | Goal | Anion Exchange Sorbent | Cation Exchange Sorbent |
|----------------|------------------------------------------------|--------------------------------------------|--------------------------------------------|
| | | pH | pH |
| WASH | To promote bonding between sorbent and analyte | > Analyte pKa and / or < Sorbent pKa | < Analyte pKa and / or > Sorbent pKa |
| ELUTION | To disrupt bonding between sorbent and analyte | < Analyte pKa and / or > Sorbent pKa | > Analyte pKa and / or < Sorbent pKa |

ENVIRO-CLEAN[®] Anion Exchange Phase

AMINOPROPYL

Organic Loading = 6.65%
Surface Area = 500 m²/g
Pore Volume = 0.77 cm³/g

Average Pore Size = 60Å
Anion Exchange = 0.28 meq/g

COLUMNS

| Tube Volume (mL) | Sorbent Amount (mg) | Units per Pack | Part Number |
|------------------|---------------------|----------------|-------------|
| 3 | 200 | 50 | EUNAX123 |
| 3 | 500 | 50 | EUNAX153 |
| 6 | 500 | 50 | EUNAX156 |
| 6 | 1000 | 30 | EUNAX1M6 |
| 10 | 100 | 50 | EUNAX11Z |
| 15 | 2000 | 20 | EUNAX12M15 |

QUATERNARY AMINE WITH CHLORIDE COUNTER ION

Organic Loading = 8.40%
Surface Area = 500 m²/g

Average Pore Size = 60Å
Anion Exchange = 0.230 meq/g

COLUMNS

| Tube Volume (mL) | Sorbent Amount (mg) | Units per Pack | Part Number |
|------------------|---------------------|----------------|-------------|
| 1 | 100 | 100 | EUQAX111 |
| 3 | 200 | 50 | EUQAX123 |
| 3 | 500 | 50 | EUQAX153 |
| 6 | 500 | 50 | EUQAX156 |
| 6 | 1000 | 30 | EUQAX1M6 |
| 10 | 100 | 50 | EUQAX11Z |
| 10 | 200 | 50 | EUQAX12Z |

QUATERNARY AMINE WITH ACETATE COUNTER ION

Organic Loading = 8.40%
Surface Area = 500 m²/g

Average Pore Size = 60Å
Anion Exchange = 0.230 meq/g

COLUMNS

| Tube Volume (mL) | Sorbent Amount (mg) | Units per Pack | Part Number |
|------------------|---------------------|----------------|-------------|
| 3 | 500 | 50 | EAQAX153 |
| 6 | 1000 | 30 | EAQAX1M6 |
| 10 | 100 | 50 | EAQAX11Z |

QUATERNARY AMINE WITH HYDROXIDE COUNTER ION

Organic Loading = 8.40%
Surface Area = 500 m²/g

Average Pore Size = 60Å
Anion Exchange = 0.230 meq/g

COLUMNS

| Tube Volume (mL) | Sorbent Amount (mg) | Units per Pack | Part Number |
|------------------|---------------------|----------------|-------------|
| 3 | 500 | 50 | EAQAX153 |

QUATERNARY AMINE WITH FORMATE COUNTER ION

Organic Loading = 8.40%
Surface Area = 500 m²/g

Average Pore Size = 60Å
Anion Exchange = 0.230 meq/g

COLUMNS

| Tube Volume (mL) | Sorbent Amount (mg) | Units per Pack | Part Number |
|------------------|---------------------|----------------|-------------|
| 3 | 500 | 50 | EFQAX153 |



POLYIMINE

Organic Loading = 8.40%
Surface Area = 500 m²/g

Average Pore Size = 60Å
Anion Exchange = 0.230 meq/g

COLUMNS

| Tube Volume (mL) | Sorbent Amount (mg) | Units per Pack | Part Number |
|------------------|---------------------|----------------|-------------|
| 3 | 500 | 50 | EUPAX153 |
| 3 | 1000 | 50 | EUPAX1M3 |
| 6 | 500 | 50 | EUPAX156 |

ENVIRO-CLEAN[®] Anion Exchange Phase

PRIMARY/SECONDARY AMINE


Organic Loading = 11.1%

Average Pore Size = 60Å

Surface Area = 500 m²/g

Anion Exchange = 1.100 meq/g

| COLUMNS | | | |
|------------------|---------------------|----------------|-------------|
| Tube Volume (mL) | Sorbent Amount (mg) | Units per Pack | Part Number |
| 1 | 50 | 100 | EUPSA1L1 |
| 1 | 100 | 100 | EUPSA111 |
| 3 | 200 | 50 | EUPSA123 |
| 3 | 500 | 50 | EUPSA153 |
| 6 | 500 | 50 | EUPSA156 |
| 6 | 1000 | 30 | EUPSA1M6 |



ENVIRO-CLEAN[®] Cation Exchange Phase

CARBOXYLIC ACID

Organic Loading = 8.75%
Surface Area = 500 m²/g

Average Pore Size = 60Å
Cation Exchange = 0.043 meq/g

COLUMNS

| Tube Volume (mL) | Sorbent Amount (mg) | Units per Pack | Part Number |
|------------------|---------------------|----------------|-------------|
| 1 | 100 | 100 | EUCCX111 |
| 3 | 200 | 50 | EUCCX123 |
| 3 | 500 | 50 | EUCCX153 |
| 6 | 500 | 50 | EUCCX156 |
| 6 | 1000 | 30 | EUCCX1M6 |

TRIACETIC ACID

Organic Loading = 7.50%
Surface Area = 500 m²/g
Pore Volume = 0.77 cm³/g

Average Pore Size = 60Å
Cation Exchange = 0.10 meq/g
Anion Exchange = 0.15 meq/g

COLUMNS

| Tube Volume (mL) | Sorbent Amount (mg) | Units per Pack | Part Number |
|------------------|---------------------|----------------|-------------|
| 1 | 100 | 100 | EUCCX111 |
| 3 | 200 | 50 | EUCCX123 |
| 3 | 500 | 50 | EUCCX153 |
| 6 | 500 | 50 | EUCCX156 |
| 6 | 1000 | 30 | EUCCX1M6 |

BENZENESULFONIC ACID

Organic Loading = 10.69%
Surface Area = 500 m²/g

Average Pore Size = 60Å
Cation Exchange = 0.320 meq/g

COLUMNS

| Tube Volume (mL) | Sorbent Amount (mg) | Units per Pack | Part Number |
|------------------|---------------------|----------------|-------------|
| 3 | 100 | 50 | EUBCX113 |
| 3 | 200 | 50 | EUBCX123 |
| 3 | 500 | 50 | EUBCX153 |
| 6 | 500 | 50 | EUBCX156 |
| 6 | 1000 | 30 | EUBCX1M6 |
| 10 | 100 | 50 | EUBCX11Z |
| 10 | 200 | 50 | EUBCX12Z |
| 10 | 500 | 50 | EUBCX15Z |
| 15 | 2000 | 20 | EUBCX12M15 |

BENZENESULFONIC ACID HIGH LOAD

Organic Loading = 16.50%
Surface Area = 500 m²/g

Average Pore Size = 60Å
Cation Exchange = 0.650 meq/g

COLUMNS

| Tube Volume (mL) | Sorbent Amount (mg) | Units per Pack | Part Number |
|------------------|---------------------|----------------|-------------|
| 1 | 100 | 100 | EUBCX1HL11 |
| 3 | 50 | 50 | EUBCX1HLL3 |
| 3 | 100 | 50 | EUBCX1HL13 |
| 6 | 500 | 50 | EUBCX1HL56 |
| 6 | 1000 | 30 | EUBCX1HLM6 |
| 10 | 200 | 50 | EUBCX1HL2Z |
| 10 | 500 | 50 | EUBCX1HL5Z |

PROPYLSULFONIC ACID

Organic Loading = 7.00%
Surface Area = 500 m²/g


Average Pore Size = 60Å
Cation Exchange = 0.180 meq/g

COLUMNS

| Tube Volume (mL) | Sorbent Amount (mg) | Units per Pack | Part Number |
|------------------|---------------------|----------------|-------------|
| 3 | 500 | 50 | EUPCX153 |
| 10 | 500 | 50 | EUPCX15Z |

ENVIRO-CLEAN[®] COPOLYMERIC EXTRACTION SORBENTS

Copolymeric sorbents are composed of a silica backbone bonded with two types of functional chains. One is either an ion exchanger or polar chain, while the other is a hydrophobic carbon chain. The copolymeric phases manufactured by UCT are produced in a way as to allow for equal parts of each functional group to attach to the silica substrate yielding reproducible bonded phases and unique copolymeric chemistries. This type of mixed-mode separation is beneficial when one is looking to extract both neutral and charged compounds and typically results in a cleaner, final extract.





| Sorbent | Category | Structure | pKa |
|-----------------------------|---------------|------------------------|----------------|
| Benzenesulfonic Acid (BCX2) | Strong Cation | C8+-Si-(CH2)2-Ph-SO3- | Always Charged |
| Benzenesulfonic Acid (BCX3) | Strong Cation | C18+-Si-(CH2)2-Ph-SO3- | Always Charged |
| Propylsulfonic Acid (PCX2) | Strong Cation | C8+-Si-(CH2)3SO3H | <1 |
| Carboxylic Acid (CCX2) | Weak Cation | C8+-Si-(CH2)2COOH | 4.8 |
| Quaternary Amine (QAX2) | Strong Anion | C8+-Si-(CH2)3N+(CH3)3 | Always Charged |
| Quaternary Amine (QAX3) | Strong Anion | C18+-Si-(CH2)3N+(CH3)3 | Always Charged |
| Aminopropyl (NAX2) | Weak Anion | C8+-Si-(CH2)3NH2 | 9.8 |
| Cyanopropyl (CNP2) | Hydrophilic | C8+-Si-(CH2)3CN | N/A |
| Cyclohexyl (CYH2) | Hydrophobic | C8+-Si-(CH2)-C6H12 | N/A |

ENVIRO-CLEAN[®] Copolymeric Exchange Phase



| Analytes | Washes | Elutions |
|----------------|--------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------|
| Cations/Anions | 1) Aqueous to remove polar matrix components. | 1) Organic to elute hydrophobically bound analytes. |
| Alkanes | 2) Methanol to disrupt residual hydrophobic interferences and to remove any remaining residual matrix. | 2) Organic with a pH that would neutralize ionically bound analytes. |
| Alkenes | | 3) Aqueous buffer with high ionic strength. |
| Aromatics | | 4) Solvent possessing a counter ion that would bond to sorbent and displace analyte of interest. |

EXTRACTION MECHANISMS OF COPOLYMERIC BONDED PHASES

A sample composed of a theoretical neutral parent drug and its charged (acidic) metabolite is applied at a pH of 6 (**Figure 1**). At this pH, most amine groups are positively charged. Since this sorbent is positively charged, compounds with positively charged cations are repelled. Depending on the pKa of the metabolite, the carboxylic acid groups may be negatively charged, allowing the metabolite to bond to the positively charged sorbent. The column also possesses a hydrophobic chain which allows the neutral parent drug to bond to the sorbent. Water or a weak aqueous buffer (pH 6) washes away hydrophilically bound interferences (**Figure 2**). The column is then dried utilizing organic solvent to ensure it's free of any residual aqueous phase that would hinder effective elution. After drying, analytes of interest can be eluted using a two-step process. During the first elution (**Figure 3**), the hydrophobically bound neutral parent drug is eluted with a solvent of minimal polarity, such as hexane/ethyl acetate (80:20). The second elution (**Figure 4**) employs an acid to neutralize the charge of acidic analytes. The ionic interaction is released, and analytes are eluted in an appropriate solvent mixture. If a tiered elution scheme is not desired, a universal solvent can be put to use that can effectively disrupt both bonding mechanisms at once to elute all analytes of interest within one single elution step.



ENVIRO-CLEAN[®] Copolymeric Exchange Phase



C8 PLUS CYCLOHEXYL

Organic Loading = 14.0%

Average Pore Size = 60Å

Surface Area = 500 m²/g

Pore Volume = 0.77 cm³/g

COLUMNS

| Tube Volume (mL) | Sorbent Amount (mg) | Units per Pack | Part Number |
|------------------|---------------------|----------------|-------------|
| 3 | 200 | 50 | EUCYH223 |
| 6 | 500 | 50 | EUCYH256 |

C8 PLUS CYANOPROPYL

Organic Loading = 14.0%

Average Pore Size = 60Å

Surface Area = 500 m²/g

Pore Volume = 0.77 cm³/g

COLUMNS

| Tube Volume (mL) | Sorbent Amount (mg) | Units per Pack | Part Number |
|------------------|---------------------|----------------|-------------|
| 3 | 200 | 50 | EUCYH223 |

C8 PLUS PROPYLSULFONIC

Organic Loading = 14.62%

Average Pore Size = 60Å

Surface Area = 500 m²/g

Pore Volume = 0.77 cm³/g

Exchange Capacity = 0.11 meq/g

COLUMNS

| Tube Volume (mL) | Sorbent Amount (mg) | Units per Pack | Part Number |
|------------------|---------------------|----------------|-------------|
| 3 | 200 | 50 | EUPCX223 |

C8 PLUS CARBOXYLIC ACID

Organic Loading = 11.45%

Average Pore Size = 60Å

Surface Area = 500 m²/g

Exchange Capacity = 0.110 meq/g

COLUMNS

| Tube Volume (mL) | Sorbent Amount (mg) | Units per Pack | Part Number |
|------------------|---------------------|----------------|-------------|
| 3 | 200 | 50 | EUCCX223 |
| 6 | 1000 | 30 | EUCCX2M6 |



ENVIRO-CLEAN[®] Copolymeric Exchange Phase

C8 PLUS BENZENESULFONIC ACID

Organic Loading = 12.40%
Surface Area = 500 m²/g
Pore Volume = 0.77 cm³/g

Average Pore Size = 60Å
Exchange Capacity = 0.077 meq/g

| COLUMNS | | | |
|------------------|---------------------|----------------|-------------|
| Tube Volume (mL) | Sorbent Amount (mg) | Units per Pack | Part Number |
| 3 | 200 | 50 | EUBCX223 |
| 6 | 500 | 50 | EUBCX256 |
| 6 | 1000 | 30 | EUBCX2M6 |
| 10 | 200 | 50 | EUBCX22Z |
| 10 | 500 | 50 | EUBCX25Z |

C8 PLUS AMINOPROPYL

Organic Loading = 12.10%
Surface Area = 500 m²/g
Pore Volume = 0.77 cm³/g

Average Pore Size = 60Å
Exchange Capacity = 0.144 meq/g

| COLUMNS | | | |
|------------------|---------------------|----------------|-------------|
| Tube Volume (mL) | Sorbent Amount (mg) | Units per Pack | Part Number |
| 3 | 200 | 50 | EUNAX223 |
| 3 | 500 | 50 | EUNAX253 |
| 6 | 1000 | 30 | EUNAX2M6 |
| 10 | 100 | 50 | EUNAX21Z |
| 75 | 1000 | 10 | EUNAX210M75 |

C18 PLUS BENZENESULFONIC ACID

Organic Loading = 12.4%
Surface Area = 500 m²/g

Average Pore Size = 60Å
Exchange Capacity = 0.077 meq/g

| COLUMNS | | | |
|------------------|---------------------|----------------|-------------|
| Tube Volume (mL) | Sorbent Amount (mg) | Units per Pack | Part Number |
| 3 | 50 | 50 | EUBCX3L3 |
| 3 | 100 | 50 | EUBCX313 |

C8 PLUS QUATERNARY AMINE

Organic Loading = 13.00%
Surface Area = 500 m²/g

Average Pore Size = 60Å
Exchange Capacity = 0.170 meq/g

| COLUMNS | | | |
|------------------|---------------------|----------------|-------------|
| Tube Volume (mL) | Sorbent Amount (mg) | Units per Pack | Part Number |
| 3 | 200 | 50 | EUQAX223 |
| 3 | 500 | 50 | EUQAX253 |
| 6 | 500 | 50 | EUQAX256 |
| 6 | 1000 | 30 | EUQAX2M6 |

ENVIRO-CLEAN[®] Polymeric Exchange Phase

ENVIRO-CLEAN[®] POLYMERIC DVB EXTRACTION SORBENT

Enviro-Clean[®] polymeric extraction sorbents are formulated with an ultra- clean, highly cross-linked styrene and divinylbenzene polymer sorbent. The sorbent can be functionalized with many of the same phases as our silica based sorbents. Possibilities include standard ion exchange functionalities. Enviro-Clean[®] particles have an average particle size of 30 microns with enhanced loading capacity. This higher capacity translates into lower requirements for bed mass. The Enviro-Clean[®] sorbent also eliminates the need for the initial column conditioning step. All of these attributes ultimately result in improved cost to the end user.

ADVANTAGES OF ENVIRO-CLEAN[®] DVB SORBENT:

- No conditioning steps
- High and reproducible recoveries
- Cross-linked sorbent minimized bead swelling
 - Reduced sorbent mass
 - Improved flow rates
 - pH stable from 1-14
 - High sorbent capacity



| COLUMNS | | | | |
|------------------|---------------------|-----------|----------------|-------------|
| Tube Volume (mL) | Sorbent Amount (mg) | Frit Type | Units per Pack | Part Number |
| 3 | 50 | PTFE | 50 | ECDVB1L3 |
| 3 | 60 | PE | 100 | ECDVB063P |
| 3 | 100 | PTFE | 50 | ECDVB113 |
| 6 | 200 | PE | 30 | ECDVB126P |
| 6 | 500 | PTFE | 30 | ECDVB156 |
| 6 | 500 | PE | 30 | ECDVB156P |

ENVIRO-CLEAN[®] Glass Cartridge


ENVIRO-CLEAN[®] INERT GLASS CARTRIDGE

For ultra-clean extractions, UCT offers inert, 6 mL glass cartridges in a variety of bonded phases packed between two PTFE frits.



| COLUMNS | | | | | |
|------------------|---------------------|--------------------------------|-----------|----------------|-------------|
| Tube Volume (mL) | Sorbent Amount (mg) | Bonded Phase | Endcapped | Units per Pack | Part Number |
| 6 | 500 | C2-ETHYL | YES | 100 | EEC02156G |
| 6 | 500 | C8-OCTYL | YES | 50 | EEC08156G |
| 6 | 500 | C18-OCTADECYL | NO | 50 | EUC18156G |
| 6 | 500 | C18-OCTADECYL | YES | 50 | EEC18156G |
| 6 | 1000 | C18-OCTADECYL | NO | 30 | EUC181M6G |
| 6 | 1000 | C18-OCTADECYL | YES | 50 | EEC181M6G |
| 6 | 1000 | CYCLOHEXL | YES | 50 | EECYH1M6G |
| 6 | 1000 | PHENYL | NO | 50 | EUPHY1M6G |
| 6 | 1000 | ALUMINA-NEUTRAL | N/A | 30 | EUALN1M6G |
| 6 | 500 | SILICA-ACID WASHED | N/A | 10 | EUSIL156G |
| 6 | 1000 | SILICA-ACID WASHED | N/A | 10 | EUSIL1M6G |
| 6 | 500 | PHARMA-SIL [®] SILICA | N/A | 10 | EPHSIL156G |
| 6 | 200 | POLYSTYRENE DVB | N/A | 10 | ECDVB126G |

ENVIRO-CLEAN® Bulk Sorbent Guide



ENVIRO-CLEAN® BULK SORBENT GUIDE



UCT sorbents are available in bulk quantities in a variety of sizes ranging from 10g to 50Kg. Common bulk sorbents used in environmental analysis are listed below, but all varieties of bonded phases can be packaged in comparable quantities.

| PRODUCT NAME | PART NUMBER | QUANTITY |
|-----------------------------------------------------------------------------------------|-------------|------------------|
| Florisil® A 100-200 Mesh | ECFLOR00D | 500g |
| | ECFLOR03K | 3kg |
| Florisil® PR | ECFLS00D | 500g |
| | ECFLS03K | 3kg |
| Silica Gel 100-200 Mesh suitable for column chromatography | ECSIOH00D | 500g |
| | ECSIOH03K | 3kg |
| Sodium Sulfate Anhydrous ACS Grade Granular 60 Mesh | ECSS01K | 1kg |
| | ECSS05K | 5kg |
| | ECSS10K | 10kg |
| | ECSS25K | 25kg |
| Magnesium Sulfate Anhydrous Organic Free Powder Reagent, 99.5% min. | ECMAG00C | 100g (1 unit) |
| | ECMAG00D | 500g (1 unit) |
| | ECMAG00DCS | 1 case (4 units) |
| Celite 566 Hydromatrix Substitute | EC56601K | 1kg |
| | EC56603K | 3kg |
| | EC566030K | 30kg |
| Ottawa Sand | ECOTT01K | 1kg |
| | ECOTT05K | 5kg |
| | ECOTT10K | 10kg |
| | ECOTT25K | 25kg |
| Graphitized Carbon 120/400 Mesh | EUCARB00C | 100g |
| | EUCARB00K | 1kg |



HydraFlow®

4 Channel Disc/Cartridge Manifold



HYDRAFLOW[®] 4 CHANNEL CARTRIDGE MANIFOLD

FEATURES

- **Liquid Channel Switching** – Patented channel-switching design to ensure the convenient and efficient transition between sample collection and waste discharge. This optimized feature improves the overall user experience and robustness of the manifold by eliminating the need for reopening the vacuum chamber mid-extraction, repeated on/off switching of the vacuum pump, time-consuming venting and relocation of collection bottles, and the cumbersome draining of the waste liquid and voiding of the vacuum chamber.
- **Multi-Sample Processing** – Four channels can be used individually or simultaneously based on the user preference.
- **No Need for Glass Cartridge Adapters** – HydraFlow[®] does not require the use of expensive glass cartridge adapters which too often break with prolonged usage and lead to costly replacement.
- **Precise Flow Control Capabilities** – Featured control valves ensure accurate sample flow rates ranging from 1 mL/min to 45 mL/min and high reproducibility across a wide range of extractions.
- **Rugged Anti-Corrosion Design** – Manifold parts feature PTFE or stainless-steel composition which can resist degradation from prolonged exposure to organic solvents, such as Dichloromethane, and acids. This provides enhanced chemical resistance and long-term durability to the unit.
- **Liquid Level Visualization** – Graduated tick marks on Universal cartridge adapters allow for staged liquid-level visualization during an extraction.
- **Separation of Organic and Aqueous Waste** – Separate channels for aqueous waste collection and organic waste collection contribute to long-term cost savings when it comes to waste disposal.
- **Lightweight and Compact Footprint** – Simple design allows for fume hood setup if desired and easy relocating when necessary, throughout the lab.

TECHNICAL SPECIFICATIONS

Dimensions: 520 mm x 290 mm x 440 mm (L x W x H)

Weight: 12 Kgs (26 lbs)

Vacuum Requirements: 25" Hg (minimum)

Collection Devices Used: 40 mL Collection Vials

Sample Capacity Per Unit: 1-4 Samples



HYDRAFLOW[®] UNIT OVERVIEW

COMPLETE CONFIGURATION



4 x Front Flow Path Valves – These valves control individual channels and are used to switch between Waste Flow, Elution Flow, and OFF position.

1 x Waste Separator Valve – Contains 3 varying positions: Aqueous, Waste (OFF) and Organic to segregate and simultaneously collect different solvent classes.

Aqueous – Used when processing water samples. This valve channels all the waste to the 20L Aqueous Waste Collection Container.

Waste (OFF) – Used to turn off the vacuum.

Organic – Used during pre-rinsing and/or pre-conditioning the SPE cartridge. This valve channels all the waste to the 1L Organic Waste Collection Bottle.

2 x Fine-Tuning Flow Control Knobs – The front solvent channel switches can be further refined using these independent controls for Waste Flow and Elution Flow.

Waste Flow Control Knob – Allows for the simultaneous flow regulation of all 4 channels when the waste is being collected during pre-conditioning, sample addition and wash procedure. When using this control knob, the Front Flow Path Valve should be set to the “Waste Flow” to begin.

Elution Flow Control Knob – Allows for the simultaneous flow regulation of all 4 channels during the elution step into the sample collection vials. When using this control knob, the Front Flow Path Valve should be set to the “Elution Flow” to begin.

MANIFOLDS



ENVIRO-CLEAN® DISK MANIFOLD

DISK MANIFOLD AND ACCESSORIES



6 Station Manifold



Glass Cartridge Adapter




Bottle Holder



Jar Holder



Universal Cartridge Adapter



47mm Aluminum Clamp




90mm Aluminum Clamp



47mm Support Base



90mm Support Base



47mm 300ml Funnel



90mm 1000ml Funnel



47mm/90mm KEL-F Screen



Vacuum Pump



Teflon Stopcock and Body



Waste Trap

| Description | Units | Part Number |
|----------------------------------------------------------------------------------|-------|------------------|
| Manifolds | | |
| 1 Station Manifold | 1 | ECUCTVAC1 |
| 3 Station Manifold | 1 | ECUCTVAC3 |
| 6 Station Manifold | 1 | ECUCTVAC6 |
| NOTE: The above numbers indicate the stainless steel manifolds base only. | | |
| UCT Universal Cartridge Accessories | | |
| Glass cartridge adapter | 1 | ECUCTADP |
| Universal cartridge bottle holder adapter | 1 | ECUNIBHD-PP |
| Universal cartridge jar bottle holder adapter | 1 | ECUNIJD100-PP |
| Universal cartridge adapter (Compatible with JT Baker Manifold*) | 1 | ECBMADP |
| SPE Disk Accessories | | |
| 47mm aluminum clamp | 1 | ECCG1420 |
| 90mm aluminum clamp | 1 | ECUC0502 |
| 47mm support base | 1 | ECQSB47 |
| 90mm support base | 1 | ECQSB90 |
| 47mm 300 mL funnel | 1 | ECQFN300 |
| 90mm 1000 mL funnel | 1 | ECQFN1000 |
| 47mm KEL-F screen | 1 | ECUCT47 |
| 90mm KEL-F screen | 1 | ECUCT90 |
| Additional Accessories | | |
| Vacuum pump - 110 volt | 1 | ECROCKER400 |
| Vacuum pump - 220 volt | 1 | ECROCKER400-220V |
| Teflon stopcock and body | 1 | ECUCTSC |
| Waste trap | 1 | ECUCTTRAP20 |

ENVIRO-CLEAN® GLASS BLOCK MANIFOLD

GLASS BLOCK MANIFOLD

A complete Vacuum Manifold System consists of a glass block, Corian® manifold lid, a cover gasket, vacuum gauge and assembly, PTFE tips, an adjustable collection rack, bulkhead luer fittings, plugs and a glass block safety tray. The Vacuum Manifold System is available in either 16 or 24 positions.

These manifold systems are durable and highly chemical resistant, designed to provide years of trouble-free extractions.



| Description | Part Number |
|---------------------------------------------|-------------|
| Complete 16 Position Vacuum Manifold System | VMF016GL |
| Complete 24 Position Vacuum Manifold System | VMF024GL |

ENVIRO-CLEAN® GLASS BLOCK MANIFOLD

GLASS BLOCK MANIFOLD ACCESSORIES




Glass Block




Manifold Lid Legs



Gasket



Collection Rack



Collection Rack
(12 Position)



Collection Rack Posts



Collection Rack
Retaining Clips



Vacuum Gauge and
Bleed Valve



Bulkhead Luer Fittings



Luer Plugs



Flange Caps



Large Volume Transfer
Tubes



Cartridge Adapters



20L Waste Trap

| Description | Units | Part Number |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------|-------------------------|
| Manifold Lid Legs – The lid legs can be used to set the manifold lid on a surface while loading columns, changing collection tubes or removing waste. | 4 | VMF02120-1 |
| Gasket – A foam gasket that fits both the 16 and 24 position lids. | 2 | VMF04121 |
| Collection Rack (16 position) – A polypropylene rack that is highly resistant to chemical degradation and abuse. This rack allows the use of 13 and 16 mm disposable test tubes. | 1 | VMF06125 |
| Collection Rack (24 position) – A polypropylene rack that is highly resistant to chemical degradation and abuse. This rack allows the use of 13 and 16 mm disposable test tubes. | 1 | VMF04125 |
| Collection Rack (12 position) – A polypropylene rack that is highly resistant to chemical degradation and abuse. This rack is designed for the use of 27 mm (VOA vials) and smaller disposable collection vials. | 1 | VMF02125 |
| Collection Rack Posts – These posts can be ordered as replacement parts for the posts in all collection racks. | 3 | VMF02127 |
| Collection Rack Retaining Clips – These clips are replacement parts for the clips included in all collection racks. | 12 | VMF02129 |
| Vacuum Gauge and Bleed Valve – This system is used in monitoring and adjusting vacuum. | 1 | VMF02122 |
| Bulkhead Luer Fittings – These fittings screw into the lid allowing the sample to transfer from the column into the PTFE Luer tip to the test tube. | 12 | VMF21BFN |
| Luer Plugs – These plugs fit into the bulkhead fittings in order to seal unused bulkhead fittings. These can also be used to break vacuum to the manifold. | 12 | VMF21PLN |
| | 1 mL | 50 |
| | 3 mL | 50 |
| Flange Caps – Used with the Luer Caps, Flange Caps plug the top of SPE cartridges. | 6 & 10 mL | 50 |
| | 15 mL | 50 |
| | 25 mL | 50 |
| | | CR0001P |
| | | CR0004P |
| | | CR0008P |
| | | CR0015P |
| | | CR0025P |
| Large Volume Transfer Tubes – Used to transfer large volumes (100-1000mL) from a water collection bottle to an SPE cartridge. | 6 | VMFSTFR06 |
| | 12 | VMFSTFR12 |
| Large Volume Transfer Tubes For Perfluorinated Compound Analysis – Used to transfer large volumes (100-1000mL) from a water collection bottle to an SPE cartridge. | 6 | |
| | 12 | |
| 1, 3, 6, 10 mL Cartridge Adapters | 15 | AD0000AS |
| | 100 | AD00000C |
| 20L Waste Trap | 1 | ECUCTTRAP20 |
| 20L Waste Trap Adapter – 3/8" x 1/4" PVDF ADPT for fitting to glass block manifold. | 1 | ECUCTTRAP20-ADPT |

SELECTRA® U/HPLC Columns



SELECTRA[®] U/HPLC COLUMNS



SELECTRA[®] DA

- Unique Biphenyl phase.
- Excellent selectivity for a wide range of compounds including therapeutic drugs, drugs of abuse, mycotoxins, veterinary drugs and pesticides.
- Ability to retain compounds that can be difficult to retain on a C18
- Can achieve significant selectivity changes with the choice of acetonitrile or methanol as the organic solvent.
- Carbon Load 13%.
- Conforms to USP L11.
- Fully endcapped.

The SELECTRA[®] line of HPLC columns is created using an ultra-high purity, Type B, spherical silica. This support material minimizes surface activity and allows for high density functional group bonding. Columns are available with either 1.8, 3, or 5 μm particle sizes

| Guard Column Holder | |
|---------------------------------------------|------------------------|
| Description | Part Number |
| HPLC Guard Cartridge Holder | SLGRDHLDLDR |
| UHPLC Direct Connect Guard Cartridge Holder | SLGRDHLDLDR-HPOPT |
| Replacement Peek Tip for Holder | SLGRDHLDLDR-TIP (2/pk) |

| SELECTRA [®] DA | | | |
|--------------------------|------------------|-------------------|------------------|
| Column Length (mm) | Column i.d. (mm) | Particle Size | Part Number |
| 50 | 2.1 | 1.8 μm | SLDA50ID21-18UM |
| 100 | 2.1 | 1.8 μm | SLDA100ID21-18UM |
| 50 | 4.6 | 1.8 μm | SLDA50ID46-18UM |
| 100 | 4.6 | 1.8 μm | SLDA100ID46-18UM |
| 50 | 2.1 | 3 μm | SLDA50ID21-3UM |
| 100 | 2.1 | 3 μm | SLDA100ID21-3UM |
| 50 | 4.6 | 3 μm | SLDA50ID46-3UM |
| 100 | 4.6 | 3 μm | SLDA100ID46-3UM |
| 150 | 4.6 | 3 μm | SLDA150ID46-3UM |
| 50 | 2.1 | 5 μm | SLDA50ID21-5UM |
| 100 | 2.1 | 5 μm | SLDA100ID21-5UM |
| 50 | 4.6 | 5 μm | SLDA50ID46-5UM |
| 100 | 4.6 | 5 μm | SLDA100ID46-5UM |
| 150 | 4.6 | 5 μm | SLDA150ID46-5UM |

| Guard Cartridge Columns (2/pack)* | | | |
|-----------------------------------|------------------|-------------------|-------------------|
| Column Length (mm) | Column i.d. (mm) | Particle Size | Part Number |
| 10 | 2.0 | 1.8 μm | SLDAGDC20-18UMOPT |
| 10 | 2.0 | 3 μm | SLDAGDC21-3UM |
| 10 | 2.0 | 5 μm | SLDAGDC21-5UM |
| 10 | 4.6 | 1.8 μm | SLDAGDC46-18UMOPT |
| 10 | 4.6 | 3 μm | SLDAGDC46-3UM |
| 10 | 4.6 | 5 μm | SLDAGDC46-5UM |

*Guard Cartridge columns must be used with a UCT guard cartridge holder.

Storage of LC Columns:

Do not allow LC analytical columns to stand uncapped for any length of time. Store an LC column in methanol or an appropriate organic solvent, capped at both ends. A dry LC column is sometimes difficult to reactivate and may not recover to full performance status.

SELECTRA® U/HPLC COLUMNS



SELECTRA® C18

- Optimum retention for traditional reverse phase analysis.
- Highest hydrophobic interactions in the Selectra® column line.
- Carbon Load 20%.
- Conforms to USP L1.
- Fully endcapped.

SELECTRA® Aqueous C18

- Polar modified C18.
- Similar non-polar retention to traditional C18.
- Enhanced selectivity and retention for difficult to retain polar analytes.
- Suitable in up to 100% aqueous mobile phases.
- Carbon Load 10%.
- Conforms to USP L1.
- Fully endcapped.

| SELECTRA® C18 | | | |
|--------------------|------------------|---------------|--------------------|
| Column Length (mm) | Column i.d. (mm) | Particle Size | Part Number |
| 50 | 2.1 | 1.8 µm | SLC-1850ID21-18UM |
| 100 | 2.1 | 1.8 µm | SLC-18100ID21-18UM |
| 50 | 4.6 | 1.8 µm | SLC-1850ID46-18UM |
| 100 | 4.6 | 1.8 µm | SLC-18100ID46-18UM |
| 50 | 2.1 | 3 µm | SLC-1850ID21-3UM |
| 100 | 2.1 | 3 µm | SLC-18100ID21-3UM |
| 50 | 4.6 | 3 µm | SLC-1850ID46-3UM |
| 100 | 4.6 | 3 µm | SLC-18100ID46-3UM |
| 150 | 4.6 | 3 µm | SLC-18150ID46-3UM |
| 50 | 2.1 | 5 µm | SLC-1850ID21-5UM |
| 100 | 2.1 | 5 µm | SLC-18100ID21-5UM |
| 50 | 4.6 | 5 µm | SLC-1850ID46-5UM |
| 100 | 4.6 | 5 µm | SLC-18100ID46-5UM |
| 150 | 4.6 | 5 µm | SLC-18150ID46-5UM |

| Guard Cartridge Columns (2/pack)* | | | |
|-----------------------------------|-----|--------|---------------------|
| 10 | 2.0 | 1.8 µm | SLC-18GDC20-18UMOPT |
| 10 | 2.0 | 3 µm | SLC-18GDC20-3UM |
| 10 | 2.0 | 5 µm | SLC-18GDC20-5UM |
| 10 | 4.6 | 1.8 µm | SLC-18GDC46-18UMOPT |
| 10 | 4.6 | 3 µm | SLC-18GDC46-3UM |
| 10 | 4.6 | 5 µm | SLC-18GDC46-5UM |

*Guard Cartridge columns must be used with a UCT guard cartridge holder.

| SELECTRA® Aqueous C18 | | | |
|-----------------------|------------------|---------------|------------------|
| Column Length (mm) | Column i.d. (mm) | Particle Size | Part Number |
| 50 | 2.1 | 1.8 µm | SLAQ50ID21-18UM |
| 100 | 2.1 | 1.8 µm | SLAQ100ID21-18UM |
| 50 | 4.6 | 1.8 µm | SLAQ50ID46-18UM |
| 100 | 4.6 | 1.8 µm | SLAQ100ID46-18UM |
| 50 | 2.1 | 3 µm | SLAQ50ID21-3UM |
| 100 | 2.1 | 3 µm | SLAQ100ID21-3UM |
| 50 | 4.6 | 3 µm | SLAQ50ID46-3UM |
| 100 | 4.6 | 3 µm | SLAQ100ID46-3UM |
| 150 | 4.6 | 3 µm | SLAQ150ID46-3UM |
| 50 | 2.1 | 5 µm | SLAQ50ID21-5UM |
| 100 | 2.1 | 5 µm | SLAQ100ID21-5UM |
| 50 | 4.6 | 5 µm | SLAQ50ID46-5UM |
| 100 | 4.6 | 5 µm | SLAQ100ID46-5UM |
| 150 | 4.6 | 5 µm | SLAQ150ID46-5UM |

| Guard Cartridge Columns (2/pack)* | | | |
|-----------------------------------|-----|--------|-------------------|
| 10 | 2.0 | 1.8 µm | SLAQGDC20-18UMOPT |
| 10 | 2.0 | 3 µm | SLAQGDC20-3UM |
| 10 | 2.0 | 5 µm | SLAQGDC20-5UM |
| 10 | 4.6 | 1.8 µm | SLAQGDC46-18UMOPT |
| 10 | 4.6 | 3 µm | SLAQGDC46-3UM |
| 10 | 4.6 | 5 µm | SLAQGDC46-5UM |

*Guard Cartridge columns must be used with a UCT guard cartridge holder.

SELECTRA® U/HPLC COLUMNS



SELECTRA® PFPP

- Can be used for reverse phase, normal phase, or HILIC separations.
- Exhibits multiple selectivity mechanisms including hydrogen bonding, dipole dipole, pi-pi overlap, hydrophilic (HILIC), and hydrophobic interactions.
- Carbon Load 11%.
- Conforms to USP L43.
- Fully endcapped.

SELECTRA® C8

- Less retentive, less hydrophobic than standard C18 column.
- Selectivity similar to C18 for non-polar compounds.
- Carbon Load 12%.
- Conforms to USP L7.
- Fully endcapped.

| SELECTRA® PFPP | | | |
|--------------------|------------------|---------------|--------------------|
| Column Length (mm) | Column i.d. (mm) | Particle Size | Part Number |
| 50 | 2.1 | 1.8 µm | SLPFPP50ID21-18UM |
| 100 | 2.1 | 1.8 µm | SLPFPP100ID21-18UM |
| 50 | 4.6 | 1.8 µm | SLPFPP50ID46-18UM |
| 100 | 4.6 | 1.8 µm | SLPFPP100ID46-18UM |
| 50 | 2.1 | 3 µm | SLPFPP50ID21-3UM |
| 100 | 2.1 | 3 µm | SLPFPP100ID21-3UM |
| 50 | 4.6 | 3 µm | SLPFPP50ID46-3UM |
| 100 | 4.6 | 3 µm | SLPFPP100ID46-3UM |
| 150 | 4.6 | 3 µm | SLPFPP150ID46-3UM |
| 50 | 2.1 | 5 µm | SLPFPP50ID21-5UM |
| 100 | 2.1 | 5 µm | SLPFPP100ID21-5UM |
| 50 | 4.6 | 5 µm | SLPFPP50ID46-5UM |
| 100 | 4.6 | 5 µm | SLPFPP100ID46-5UM |
| 150 | 4.6 | 5 µm | SLPFPP150ID46-5UM |
| 250 | 4.6 | 5 µm | SLPFPP250ID46-5UM |

| Guard Cartridge Columns (2/pack)* | | | |
|-----------------------------------|-----|--------|---------------------|
| 10 | 2.0 | 1.8 µm | SLPFPPGDC20-18UMOPT |
| 10 | 2.0 | 3 µm | SLPFPPGDC20-3UM |
| 10 | 2.0 | 5 µm | SLPFPPGDC20-5UM |
| 10 | 4.6 | 1.8 µm | SLPFPPGDC46-18UMOPT |
| 10 | 4.6 | 3 µm | SLPFPPGDC46-3UM |
| 10 | 4.6 | 5 µm | SLPFPPGDC46-5UM |

*Guard Cartridge columns must be used with a UCT guard cartridge holder.

| SELECTRA® C8 | | | |
|--------------------|------------------|---------------|-------------------|
| Column Length (mm) | Column i.d. (mm) | Particle Size | Part Number |
| 50 | 2.1 | 1.8 µm | SLC-850ID21-18UM |
| 100 | 2.1 | 1.8 µm | SLC-8100ID21-18UM |
| 50 | 4.6 | 1.8 µm | SLC-850ID46-18UM |
| 100 | 4.6 | 1.8 µm | SLC-8100ID46-18UM |
| 50 | 2.1 | 3 µm | SLC-850ID21-3UM |
| 100 | 2.1 | 3 µm | SLC-8100ID21-3UM |
| 50 | 4.6 | 3 µm | SLC-850ID46-3UM |
| 100 | 4.6 | 3 µm | SLC-8100ID46-3UM |
| 150 | 4.6 | 3 µm | SLC-8150ID46-3UM |
| 50 | 2.1 | 5 µm | SLC-850ID21-5UM |
| 100 | 2.1 | 5 µm | SLC-8100ID46-5UM |
| 50 | 4.6 | 5 µm | SLC-850ID46-5UM |
| 100 | 4.6 | 5 µm | SLC-8100ID46-5UM |
| 150 | 4.6 | 5 µm | SLC-8150ID46-5UM |
| 250 | 4.6 | 5 µm | SLC-8250ID46-5UM |

| Guard Cartridge Columns (2/pack)* | | | |
|-----------------------------------|-----|--------|--------------------|
| 10 | 2.1 | 1.8 µm | SLC-8GDC21-18UMOPT |
| 10 | 2.1 | 3 µm | SLC-8GDC21-3UM |
| 10 | 2.1 | 5 µm | SLC-8GDC21-5UM |
| 10 | 4.6 | 1.8 µm | SLC-8GDC46-18UMOPT |
| 10 | 4.6 | 3 µm | SLC-8GDC46-18UMOPT |
| 10 | 4.6 | 5 µm | SLC-8GDC46-5UM |


*Guard Cartridge columns must be used with a UCT guard cartridge holder.



SELECTRACORE®

Core-Shell Columns

C18 | DA | PFPP



SELECTRACORE® CORE-SHELL COLUMNS



SelectraCore® C18 Core-Shell Columns

| Length | i.d. | Particle Size | Pore Size | Part Number |
|--------|--------|---------------|-----------|---------------|
| 50 mm | 2.1 mm | 2.7 µm | 90 Å | SCS27-C18521 |
| 100 mm | 2.1 mm | 2.7 µm | 90 Å | SCS27-C181021 |
| 150 mm | 2.1 mm | 2.7 µm | 90 Å | SCS27-C181521 |
| 50 mm | 4.6 mm | 2.7 µm | 90 Å | SCS27-C18546 |
| 100 mm | 4.6 mm | 2.7 µm | 90 Å | SCS27-C181046 |
| 150 mm | 4.6 mm | 2.7 µm | 90 Å | SCS27-C181546 |

SelectraCore® C18 Guard Columns (3/pk)

| Length | i.d. | Particle Size | Pore Size | Part Number |
|--------|--------|---------------|-----------|----------------|
| 5 mm | 2.1 mm | 2.7 µm | 90 Å | SCS27-C18GDC21 |
| 5 mm | 4.6 mm | 2.7 µm | 90 Å | SCS27-C18GDC46 |

SelectraCore® PFPP Core-Shell Columns

| Length | i.d. | Particle Size | Pore Size | Part Number |
|--------|--------|---------------|-----------|---------------|
| 50 mm | 2.1 mm | 2.7 µm | 90 Å | SCS27-PFP521 |
| 100 mm | 2.1 mm | 2.7 µm | 90 Å | SCS27-PFP1021 |
| 150 mm | 2.1 mm | 2.7 µm | 90 Å | SCS27-PFP1521 |
| 50 mm | 4.6 mm | 2.7 µm | 90 Å | SCS27-PFP546 |
| 100 mm | 4.6 mm | 2.7 µm | 90 Å | SCS27-PFP1046 |
| 150 mm | 4.6 mm | 2.7 µm | 90 Å | SCS27-PFP1546 |

SelectraCore® PFPP Guard Columns (3/pk)

| Length | i.d. | Particle Size | Pore Size | Part Number |
|--------|--------|---------------|-----------|----------------|
| 5 mm | 2.1 mm | 2.7 µm | 90 Å | SCS27-PFPGDC21 |
| 5 mm | 4.6 mm | 2.7 µm | 90 Å | SCS27-PFPGDC46 |

SelectraCore® DA Core-Shell Columns

| Length | i.d. | Particle Size | Pore Size | Part Number |
|--------|--------|---------------|-----------|--------------|
| 50 mm | 2.1 mm | 2.7 µm | 90 Å | SCS27-DA521 |
| 100 mm | 2.1 mm | 2.7 µm | 90 Å | SCS27-DA1021 |
| 150 mm | 2.1 mm | 2.7 µm | 90 Å | SCS27-DA1521 |
| 50 mm | 4.6 mm | 2.7 µm | 90 Å | SCS27-DA546 |
| 100 mm | 4.6 mm | 2.7 µm | 90 Å | SCS27-DA1046 |
| 150 mm | 4.6 mm | 2.7 µm | 90 Å | SCS27-DA1546 |

SelectraCore® DA Guard Columns (3/pk)

| Length | i.d. | Particle Size | Pore Size | Part Number |
|--------|--------|---------------|-----------|---------------|
| 5 mm | 2.1 mm | 2.7 µm | 90 Å | SCS27-DAGDC21 |
| 5 mm | 4.6 mm | 2.7 µm | 90 Å | SCS27-DAGDC46 |

SelectraCore® C18

- Popular Choice for Method Development
- Excellent reversed-phase retention capacity
- Applicable to a wide range of analyte polarities
- **Target** - Suitable for most acidic, basic, & neutral compounds
- **Applications** - Food, enviro, forensic, & pharmaceutical analysis

SelectraCore® DA

- Unique polyaromatic stationary phase with alternate C18 selectivity
- Reversed-phase retention based on pi-pi interaction
- Allows analysis of compounds that are difficult to resolve on C18
- **Target** - Aromatic, conjugated or unsaturated analytes
- **Applications** - Drugs of abuse, pain management drugs, pesticides & mycotoxins


SelectraCore® PFPP

- Fluorinated stationary phase with electron-deficient phenyl rings
- Multiple retention mechanisms include dipole-dipole, pi-pi & ion exchange interactions
- Can also be used for Normal phase and HILIC separations
- Alternate selectivity to C8 & C18
- **Target** - Suitable for acidic, basic & neutral compounds
- **Applications** - Panels consisting of Beta-blockers, Benzodiazepines, TCAs & Catecholamines



SPeVAP

32 and 48 Position Multi-function Solven Evaporator



SPeVAP

UCT

SPeVAP SPECIFICATIONS

ACCESSORIES

- Required space: 23" x 18"
- Suitable electrical supply
- Dedicated 20 Amp circuit / per every three SPeVAP® modules
- Purified compressed air or N2 source
- Flathead Screwdriver and Allen Wrench

CONTENTS

- Exhaust Line and Clamp
- Drain hose and Air line
- Manual and Installation Specs
- Fuse and Power Cord




TECHNICAL SPECIFICATIONS

| | |
|-------------------|-------------------------------------------------------------------------------------------------|
| Supply Voltage | 100-240 VAC. 50,60HZ |
| Power | Dedicated 20AMP circuit per 3 SPeVAP® units |
| Max Output Power | 500W / Max Amperage – 6A |
| Fuse | 5mm x 20mm, 10A Slow Blow |
| Gas Supply | Moisture Free inert Gas or compressed air filtered to 5µm < 60psi |
| Pressure | 0-60psi / Max Supply Flow – 96L/min. |
| Water Temperature | Ambient up to 90°C |
| Water Bath Volume | 6.5L / Use DI or Distilled Water |
| Interface | 7.7' touch screen |
| Bench Space | 23" x 18" |
| Ventilation | Exhaust hose that can be routed to a fume hood or other ventilation shaft |
| Paint | Acrylic Polyurethane which is resistant to Acetone, Ethyl Acetate, Cyclohexane, Hexane, and DCM |

| Part Number | Description |
|-------------------|---------------------------------------------------|
| VMFSPEVAP-32 | SPeVAP® 32 Position |
| VMFSPEVAP-48 | SPeVAP® 48 Position |
| VMFSPEVAP-2102 | Replacement Nozzle, PTFE Coated |
| VMFSPEVAP-2103 | Nozzle O-ring, Pack of 50 |
| VMFSPEVAP-TB309 | Exhaust hose, extension (multiples of 5 ft) |
| VMFSPEVAP-PN0199 | Air Filter |
| VMFSPEVAP-SC1277 | Exhaust line clamp |
| VMFSPEVAP-TB311 | Air line (multiples of 25 ft) |
| VMFSPEVAP-FT512 | Air line quick-disconnect fitting |
| VMFSPEVAP-TB310 | Drain hose |
| VMFSPEVAP-FT256 | Replacement barb fitting for water bath drain |
| VMFSPEVAP-TB312-4 | Exhaust hose coupler |
| VMFSPEVAP-FO4811 | Gasket, 48-position, 12-13mm or Autosampler Vials |
| VMFSPEVAP-F04812 | Gasket, 48-position, 15-18mm |
| VMFSPEVAP-F03211 | Gasket, 32-position, 12-13mm or Autosampler Vials |
| VMFSPEVAP-F03212 | Gasket, 32-position, 15-18mm |

ACCESSORIES | GC LINERS | RESERVOIRS | FRITS




GC LINERS



Gas Chromatograph Glass Liners manufactured by UCT are deactivated using a proprietary silane. The silane is manufactured by UCT Specialties, LLC, a leader in high purity, specialty silanes for the chromatographic industry.

| DESCRIPTION | INNER DIAMETER (mm) | OUTER DIAMETER (mm) | LENGTH (mm) | INSTRUMENT | UNITS | UCT Part Number |
|---------------------------------------------------------------------|---------------------|---------------------|-------------|-------------------|-------|----------------------|
| 2 mm Straight Split/Splitless | 2.0 | 6.5 | 78.5 | Agilent | 1 | GCL2MM |
| | | | | | 5 | GCL2MM-5 |
| | | | | | 25 | GCL2MM-25 |
| 2 mm Straight Split/Splitless with Deactivated Glass Wool | 2.0 | 6.5 | 78.5 | Agilent | 1 | GCL2MMGW |
| | | | | | 5 | GCL2MMGW-5 |
| | | | | | 25 | GCL2MMGW-25 |
| 2 mm Gooseneck Split/Splitless | 2.0 | 6.5 | 78.5 | Agilent | 1 | GCLGN2MM |
| | | | | | 5 | GCLGN2MM-5 |
| | | | | | 25 | GCLGN2MM-25 |
| 2 mm Gooseneck Split/Splitless with Deactivated Glass Wool | 2.0 | 6.5 | 78.5 | Agilent | 1 | GCLGN2MMGW |
| | | | | | 5 | GCLGN2MMGW-5 |
| | | | | | 25 | GCLGN2MMGW-25 |
| 4 mm Straight Split/Splitless | 4.0 | 6.5 | 78.5 | Agilent | 1 | GCL4MM |
| | | | | | 5 | GCL4MM-5 |
| | | | | | 25 | GCL4MM-25 |
| 4 mm Straight Split/Splitless with Deactivated Glass Wool | 4.0 | 6.5 | 78.5 | Agilent | 1 | GCL4MMGW |
| | | | | | 5 | GCL4MMGW-5 |
| | | | | | 25 | GCL4MMGW-25 |
| 4 mm Recessed Gooseneck Split/Splitless | 4.0 | 6.5 | 78.5 | Agilent | 1 | GCLRG4MM |
| | | | | | 5 | GCLRG4MM-5 |
| | | | | | 25 | GCLRG4MM-25 |
| 4 mm Recessed Gooseneck Split/Splitless with Deactivated Glass Wool | 4.0 | 6.5 | 78.5 | Agilent | 1 | GCLRG4MMGW |
| | | | | | 5 | GCLRG4MMGW-5 |
| | | | | | 25 | GCLRG4MMGW-25 |
| 4 mm Gooseneck Split/Splitless | 4.0 | 6.5 | 78.5 | Agilent | 1 | GCLGN4MM |
| | | | | | 5 | GCLGN4MM-5 |
| | | | | | 25 | GCLGN4MM-25 |
| 4 mm Gooseneck Split/Splitless with Deactivated Glass Wool | 4.0 | 6.5 | 78.5 | Agilent | 1 | GCLGN4MMGW |
| | | | | | 5 | GCLGN4MMGW-5 |
| | | | | | 25 | GCLGN4MMGW-25 |
| 3.4 mm Straight Split 1078/1079 Inlet | 3.4 | 5.0 | 54 | Varian/ Bruker | 1 | GCL3.4MM |
| | | | | | 5 | GCL3.4MM-5 |
| | | | | | 25 | GCL3.4MM-25 |
| 3.4 mm Straight Split with Frit Inserted 1078/1079 Inlet | 3.4 | 5.0 | 54 | Varian/ Bruker | 1 | GCL3.4MMFR |
| | | | | | 5 | GCL3.4MMFR-5 |
| | | | | | 25 | GCL3.4MMFR-25 |


RESERVOIRS



| Polypropylene Reservoirs | | | | |
|--------------------------|----------------|--------------|------------------------|-------------|
| Volume Capacity | Units per Pack | No. of Frits | Porosity of Frits (µm) | Part Number |
| 1 mL | 50 | 0 | N/A | RFV0001P |
| 1 mL | 50 | 1 | 10 | RFV01F1P |
| 1 mL | 50 | 2 | 10 | RFV02F1P |
| 1 mL | 50 | 1 | 20 | RFT01F1P |
| 1 mL | 50 | 2 | 20 | RFT02F1P |
| 4 mL | 50 | 0 | N/A | RFV0004P |
| 4 mL | 50 | 1 | 10 | RFV01F4P |
| 4 mL | 50 | 2 | 10 | RFV02F4P |
| 4 mL | 50 | 1 | 20 | RFT01F4P |
| 4 mL | 50 | 2 | 20 | RFT02F4P |
| 8 mL | 50 | 0 | N/A | RFV0008P |
| 8 mL | 50 | 1 | 10 | RFV01F8P |
| 8 mL | 50 | 2 | 10 | RFV02F8P |
| 8 mL | 50 | 1 | 20 | RFT01F8P |
| 8 mL | 50 | 2 | 20 | RFT02F8P |
| 10 mL | 50 | 0 | N/A | RFV0010P |
| 10 mL | 50 | 1 | 10 | RFV1F10P |
| 10 mL | 50 | 2 | 10 | RFV2F10P |
| 10 mL | 50 | 1 | 20 | RFT1F10P |
| 10 mL | 50 | 2 | 20 | RFT2F10P |
| 15 mL | 50 | 0 | N/A | RFV0015P |
| 15 mL | 50 | 1 | 10 | RFV1F15P |
| 15 mL | 50 | 2 | 10 | RFV2F15P |
| 15 mL | 50 | 1 | 20 | RFT1F15P |
| 15 mL | 50 | 2 | 20 | RFT2F15P |
| 25 mL | 50 | 0 | N/A | RFV0025P |
| 25 mL | 50 | 1 | 10 | RFV1F25P |
| 25 mL | 50 | 2 | 10 | RFV2F25P |
| 25 mL | 50 | 1 | 20 | RFT1F25P |
| 25 mL | 50 | 2 | 20 | RFT2F25P |
| 75 mL | 50 | 0 | N/A | RFV0075P |
| 75 mL | 50 | 1 | 10 | RFV1F75P |
| 75 mL | 50 | 2 | 10 | RFV2F75P |
| 75 mL | 50 | 1 | 20 | RFT1F75P |
| 75 mL | 50 | 2 | 20 | RFT2F75P |
| 150 mL | 10 | 0 | N/A | RFV00150P |
| 150 mL | 10 | 1 | 20 | RFT1F150P |
| 150 mL | 10 | 2 | 20 | RFT2F150P |

| Glass Reservoirs | | | | |
|------------------|----------------|--------------|------------------------|-------------|
| Volume Capacity | Units per Pack | No. of Frits | Porosity of Frits (µm) | Part Number |
| 8 mL | 30 | 0 | N/A | RFV0008G |
| 8 mL | 30 | 1 | 10 | RFV01F8G |

FRITS



Polyethylene Frits

| Column Size | Diameter | Porosity | Thickness | Units | Part Number |
|-------------|----------|----------|-----------|-------|-------------|
| 1 mL | 0.232" | 10 µm | 1/16" | 100 | FR10011P |
| 1 mL | 0.232" | 20 µm | 1/16" | 100 | FR20011P |
| 1 mL | 0.232" | 20 µm | 1/8" | 100 | FT20011P |
| 4 mL | 0.357" | 7 µm | 1/16" | 100 | FR07041P |
| 4 mL | 0.357" | 10 µm | 1/16" | 100 | FR10041P |
| 4 mL | 0.357" | 20 µm | 1/16" | 100 | FR20041P |
| 4 mL | 0.357" | 20 µm | 1/8" | 100 | FT20041P |
| 4 mL | 0.357" | 100 µm | 1/16" | 100 | FR100041P |
| 8 mL | 0.498" | 10 µm | 1/16" | 100 | FR10081P |
| 8 mL | 0.498" | 20 µm | 1/16" | 100 | FR20081P |
| 8 mL | 0.513" | 20 µm | 1/8" | 100 | FT20081P |
| 10 mL | 0.357" | 10 µm | 1/16" | 100 | FR10101P |
| 10 mL | 0.357" | 20 µm | 1/16" | 100 | FR20101P |
| 10 mL | 0.357" | 20 µm | 1/8" | 100 | FT20101P |
| 15 mL | 0.630" | 10 µm | 1/16" | 100 | FR10151P |
| 15 mL | 0.641" | 20 µm | 1/16" | 100 | FR20151P |
| 15 mL | 0.641" | 20 µm | 1/8" | 100 | FT20151P |
| 25 mL | 0.792" | 10 µm | 1/16" | 100 | FR10251P |
| 25 mL | 0.816" | 20 µm | 1/8" | 100 | FT20251P |
| 75 mL | 1.050" | 10 µm | 1/16" | 100 | FR10751P |
| 75 mL | 1.050" | 20 µm | 1/8" | 100 | FT20751P |
| 150 mL | 1.515" | 20 µm | 1/8" | 20 | FT201501P |

PTFE Frits

| FR10081T | Diameter | Porosity | Thickness | Units | Part Number |
|----------|----------|----------|-----------|-------|-------------|
| 4 mL | 0.357" | 10 µm | 1.5 mm | 60 | FR10041T |
| 8 mL | 0.498" | 10 µm | 1.5 mm | 60 | FR10081T |
| 8 mL | 0.498" | 50 µm | 1.5 mm | 60 | FR10081T |
| 15 mL | 0.641" | 10 µm | 1.5 mm | 60 | FR10151T |
| 15 mL | 0.641" | 50 µm | 1.5 mm | 100 | FR50151T |

CUSTOMER SERVICE

PRICES AND TERMS

Our prices are subject to change without notice. The price in effect when we receive your order will apply. All prices are in US Dollars and are F.O.B. Lewistown, PA 17044. Terms of payment are net 30 days.

MINIMUM ORDERS

We welcome all orders, therefore, we do not have a minimum order requirement. When ordering, please include your purchase order number, complete "Ship To" and "Bill To" address, catalog number, quantity, and description of product(s). Also include your name and a phone number where you can be reached should we have any questions concerning your order.

SHIPMENTS

Normal processing is within 24 hours after receipt of an order. Unless special shipping requests have been made, our trained staff will send all orders by UPS Ground service. The appropriate shipping charges (freight & insurance costs) will be added to the invoice, unless otherwise instructed by the customer.

SPECIAL PRICING

We offer special pricing for volume purchases and standing orders. These discounts apply to bonded phase extraction column purchases only. Please call a sales representative for more information on special pricing qualifications.

RETURN POLICY

Our Quality Manager will handle all returns. Before returning merchandise, please call to obtain a return authorization number from the quality manager. We will need to know the reason for the return, date of purchase, purchase order number and invoice number in order to issue a return authorization number. Return merchandise must be received before a credit can be issued. Returns will not be accepted after 90 days. A restocking fee of 25% of the price paid, or a minimum of \$25.00 (whichever is greater) will be charged on all returns.

WARRANTY

All products manufactured by UCT are guaranteed against defects in materials and workmanship for a period of 90 days after shipment. UCT will replace any items that prove to be defective during this time period. The exclusive remedy requires the end user to first advise UCT of the defective product by phone or in writing and must include order number, the lot number and the shipping date.

To initiate this action, photographs of the product, including packaging and labeling of the containers, must be submitted to the UCT Representative for approval. With approval a return authorization can be initiated, and must be received within 30 days. Once the materials arrive at UCT a further inspection of the materials must be completed and accepted by our Quality Manager prior to further action of credits or replacement. UCT's total liability is limited to the replacement cost of UCT products.

This warranty does not apply to damage resulting from misuse.

Placing An Order

Email: info@unitedchem.com

Web: www.unitedchem.com



UCT, Inc.
2731 Bartram Rd.
Bristol, PA 19007

Phone: 800.385.3153
Fax: 215.785.1226
www.unitedchem.com



